
seca 676

	D	Bedienungsanleitung und Garantieerklärung 3
(B	Instruction manual and guarantee 50
(R	Mode d'emploi et garantie96
	IT	Manuale di istruzioni e garanzia
(S	Manual de instrucciones y garantia 190
(T	Instruções de utilização e declaração de garantia 236
(I	PL)	Instrukcja obsługi i gwarancja 282

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TABLE OF CONTENTS

1. Signed and sealed 52	6. Operation 70
2. Device description 53	6.1 Making the scale operational . 70
2.1 Congratulations! 53 2.2 Intended use 53 2.3 Description of function 53 2.4 User qualification 54 Assembly 54 Operation 54	Moving the scale
3. Safety information 55	Determining extra weight
3.1 Safety information in these instructions for use	(TARE)
4. Overview	6.3 Other functions (menu) 76 Navigating in the menu 76
4.1 Controls	Deleting saved values automatically (ACIr)
5. Before you really get started 66	Activating beeps (BEEP) 79
5.1 Scope of delivery	Setting filtering (Fil)


English

7. The seca 360° wireless network 82	8. Hygiene treatment 88
7.1 Introduction	8.1 Cleaning 88 8.2 Disinfecting 88 8.3 Sterilizing 88
Device detection 83	9. Function check 89
7.2 Operating the scale in a wireless group (menu) 84	10. Troubleshooting 89
Activating the wireless	11. Servicing
module (SYS) 84 Setting up a wireless	12. Technical data 92
group (Lrn) 84	13. Optional accessories 93
Activating automatic send	14. Spare parts 93
(ASend)	15. Disposal 94
Set the time (Time) 87	16. Warranty 94

1. SIGNED AND SEALED

With seca products, you are not only purchasing technology which has been perfected over 100 years, but also quality with official legal certification confirmed by institutions. seca products comply with European directives, standards and national laws. When you buy seca, you buy the future.

Products bearing this symbol meet the applicable regulatory requirements of the European Community, especially the following:

• Directive 93/42/EEC for medical devices

seca's professionalism is also recognized by official testing agencies. TÜV Süd Product Service, the appointed office for medical devices, confirms with this certificate that as a medical devices manufacturer, seca consistently complies with the strict legal requirements. seca's quality assurance system includes the areas of design, development, production, sales and service of medical scales and length measuring systems as well as software and measuring systems for assessing state of health and nutrition.

seca helps the environment. Saving natural resources is very important to us. We therefore make every effort to save on packaging materials wherever it makes sense and whatever is left over can be conveniently disposed of on site via the recycling system.

2. DEVICE DESCRIPTION

2.1 Congratulations!

With the **seca 676** electronic platform scale, you have purchased a highly precise and simultaneously robust device.

For over 170 years, seca has used its experience in the service of healthcare and, as a market leader, it has always set standards in many countries of the world with innovative developments in weighing and measuring.

2.2 Intended use

Electronic platform scale **seca 676** is mainly used in hospitals, medical practices and in-patient care facilities in accordance with national regulations.

The scale is for conventional determination of weight and establishment of general state of nutrition; it assists the physician supervising treatment in making a diagnosis or deciding on a course of treatment.

To make an accurate diagnosis, however, other specific examinations have to be ordered by the physician and their results taken into account, in addition to determining a weight value.

2.3 Description of function

On some variants of this scale, the weight display can be switched between kilograms (kg), pounds (lbs) and stones (sts). Weight is determined within a few seconds.

In addition to conventional determination of weight, the **seca 676** provides a function for determining body mass index. To this end, height is entered on the keypad and the body mass index associated with the weight value is calculated automatically. Length measuring devices from the **seca 360° wireless** system can transmit height to the **seca 676** wirelessly.

Via the **seca 360° wireless** wireless network, measuring results can be transmitted wirelessly to a seca wireless printer or to a PC equipped with **seca analytics** PC software and the seca USB wireless adapter.

The **seca 676** can be moved on casters.

Use the scale only for the purpose named in the section entitled "Intended use" on page 53.

2.4 User qualification

Assembly

Devices that are shipped partially assembled may only be mounted by sufficiently qualified persons such as specialist dealers, hospital technicians or seca service technicians.

Operation

The device may only be operated by healthcare professionals.

3. SAFETY INFORMATION

3.1 Safety information in these instructions for use

DANGER!

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries will occur.

WARNING!

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries may result.

CAUTION!

Used to identify a hazardous situation. If you fail to take note of this information, minor to moderate injuries may result.

ATTENTION!

Used to identify possible incorrect usage of the device. If you fail to take note of this information, you may damage the device, or the measured results may be incorrect.

NOTE:

Includes additional information about use of the device.

3.2 Basic safety instructions

Handling the instrument

- Please take note of the information in these instructions for use.
- Keep the instructions for use in a safe place. The instructions for use are a component of the device and must be available at all times.

DANGER!

Risk of explosion

Do not use the device in an environment in which one of the following gases has accumulated:

- oxygen
- flammable anesthetics
- other flammable substances/air mixtures

Λ

CAUTION!

Patient hazard, damage to device

- Additional devices which are connected to electrical medical devices must provide evidence of compliance with the relevant IEC or ISO standards (e.g. IEC 60950 for dataprocessing devices). Furthermore, all configurations must comply with the requirements of standards for medical systems (see IEC 60601-1-1 or Section 16 of the 3rd edition of IEC 60601-1 respectively). Anyone connecting additional devices to electrical medical devices is considered a system configurer and is therefore responsible for ensuring that the system complies with the requirements of standards for systems. Your attention is drawn to the fact that local laws take precedence over the above-mentioned requirements of standards. In the event of any queries, please contact your local specialist dealer or Technical Service.
- Have servicing carried out regularly as described in the relevant section of this document.
- Technical modifications may not be made to the device. The device does not contain any parts for servicing by the user. Only have servicing and repairs performed by an authorized seca service partner. You can find service partners in your area at www.seca.com or by sending an e-mail to service@seca.com.
- Only use original seca accessories and spare parts, otherwise seca will not grant any warranty.

CAUTION!

Patient hazard, malfunction

- Keep other electrical medical devices, e.g. high-frequency surgical devices, a minimum distance of approx. 1 meter away to prevent incorrect measurements or wireless transmission interference.
- Keep HF devices such as cell phones a minimum distance of approx. 1 meter away to prevent incorrect measurements or wireless transmission interference.

 The actual transmission output of HF equipment may require minimum distances of more than 1 meter. Details can be found at www.seca.com.

Preventing electric shock

WARNING!

Electric shock

- Set up the device so that the power supply socket is easy to reach and the device can be disconnected from the power supply quickly.
- Ensure that your local power supply matches the information on the power supply unit.
- Do not touch the power supply unit with wet hands.
- Do not use extension cables or power strips.
- Ensure that the power cable is not pinched or damaged by a sharp edge.
- Do not operate the device above an altitude of 3000 m.

Avoiding infections

WARNING!

Risk of infection

- Hygienically reprocess the scales regularly as described in the respective section in this document.
- Make sure that the patient has no infectious diseases.
- Make sure that the patient has no open wounds or infectious skin alterations, which may come into contact with the device.

Avoiding injuries

WARNING!

Risk of falling

- Ensure that the device is positioned firmly and level.
- Route the mains cable so that there is no risk of tripping for the patient or user.
- The device is not designed as a standing aid. Assist people with limited motor skills when they are getting up, e.g. from a wheelchair.
- Make sure that the patient does not step onto the weighing platform directly at the edges.

 Make sure that the patient steps onto the weighing platform slowly and securely.

WARNING!

Danger of slipping

- Ensure that the weighing platform is dry before the patient steps onto it.
- Ensure that the patients feet are dry before he or she steps onto the weighing platform.
- Make sure that the patient steps onto the weighing platform slowly and securely.

Preventing device damage

ATTENTION!

Damage to device

- Ensure that no liquids enter the device. They can damage the electronics.
- Switch off the device before disconnecting the power supply unit from the mains socket.
- Disconnect the power supply unit from the mains socket if you intend to not use the device for a longer period of time. Only this way it can be ensured that the device is currentless
- Make sure not to drop the device.
- Do not expose the device to any impacts or vibrations.
- Perform function controls regularly as described in the relevant section in this document. Do not operate the device if it is damaged or not working properly.
- Ensure that there is no heat source in the immediate vicinity. Do not expose to direct sunlight. The excessive temperature could damage the electronics.
- Avoid rapid temperature fluctuations. When the device is transported so that a temperature difference of more than 20 °C occurs, it must stay turned off for at least 2 hours before it can be turned on again. Otherwise, condensation water will form which can damage the electronics.
- Caustic detergents may damage the surfaces.
 Only use a soft cloth dampened with mild soapsuds to clean the surfaces of the device.

Handling measuring results

WARNING!

This device is **no** diagnostic device. It simply assists the treating physician in establishing a diagnosis.

- In order to make a precise diagnosis and initiate therapeutic measures, besides determination of the weight, further targeted examinations must be set up by the physician, and their results must be considered.
- The responsibility for diagnosis and treatment lies with the treating physician.

CAUTION!

Patient hazard

In order to avoid misinterpretations, test results for medical use must be displayed and used in SI units only. Some devices offer the ability to display test results in other units. This is only an additional function.

- Use the results exclusively in SI units.
- The use of measurement results in non-SI units is the sole responsibility of the user.

ATTENTION!

Loss of data

- Before you save measurements aguired with this device and process them for further use (e.g. in the PC software **seca analytics** or in a hospital information system), make sure the measurement values are plausible.
- If measurement values are transmitted to the seca analytics PC software or a hospital information system, make sure prior to further use that the measurement values are plausible and are assigned to the correct patient.

Handling packaging material

WARNING!

Danger of suffocation

Packaging material made of plastic foil (bags) is a choking hazard.


- Keep packaging material out of reach of children.
- In the event that the original packing material may not be available anymore, only use plastic bags with security holes in order to reduce the risk of suffocation.

NOTE:

Keep the original packing material for future use (e.g. returning for maintenance service).

4. OVERVIEW

4.1 Controls

No.	Control	Function
1	Display housing	Central control and display element, can be swiveled 180° to left and right
2	G	Switch scale on and off
3	A hold tare	Arrow key During weighing: press briefly: activate hold function press and hold: activate tare function In the menu: select submenu, select menu item increase value

No.	Control	Function
4	▼ brni menu	Arrow key During weighing: press briefly: activate BMI function press and hold: call up menu In the menu: select submenu, select menu item reduce value
5	end tring	 Enter key During weighing (if wireless network is set up): press briefly: send measuring result to devices ready to receive it (PC with USB wireless adapter) press and hold: print out measuring result (wireless printer) In the menu: confirm selected menu item save set value
6	Display	Display element for measuring results and to configure the device
7	Transport casters	The scale can be moved on these casters.
8	Locking screws	For securing the railing in the upright position
9	Access ramps	Can be folded upFacilitate access to the scale by wheelchair
10	Spirit level	Shows whether the device is horizontal
11	Power supply connection	Serves to connect the power pack supplied
12	Latch	Secures the railing in the down position during transport.
13	Foot screw	4 pcs, for precise alignment

4.2 Symbols in the display

	Symbol	Meaning
Α	\forall	Operation with power pack
В	\triangle	Non-calibratable function active
С	(1) (3) (2)	Memory location currently in use

4.3 Marking on the device and on the rating plate


Text/symbol	Meaning
Model	Model number
Ser. no.	Serial number
	Follow instructions for use
†	Electrical medical device, type B
	Insulated device, protection class II
d	Value which gives the difference between two consecutive display values
(((0123)	Device conforms to EC directives • 0123: appointed office for medical devices: TÜV Süd Product Service
F©	Symbol of the US Federal Communications Commission (FCC)
FCC ID	Device license number from the Federal Communications Commission (FCC)
IC	Device license number from Industry Canada
x-y V == max. xxx mA use compatible seca adapter only	Rating plate on the power supply connection socket • x-y V: required supply voltage • max xxx mA: maximum current consumption •
X	Do not dispose of device with household waste

4.4 Marking on the packaging

*	Protect from moisture
<u>11</u>	Arrows indicate top of product. Transport and store in an upright position.
	Fragile Do not throw or drop.
1	Permitted min. and max. temperature for transport and storage
ZZ.	Permitted min. and max. moisture for transport and storage
or crows with	Packaging material can be disposed of through recycling programs.


4.5 Menu structure

Other functions are available to you in the menu for the device. This enables you to configure the device perfectly to suit your needs (details from Page 76).

5. BEFORE YOU REALLY GET STARTED ...


5.1 Scope of delivery

No.	Component	Pcs.
Α	Scale with railing in down position	1
В	Display housing with display cable	1
С	Pan-head tapping screw B 3.5 x 9.5	3
D	Screwdriver	1
E	Power pack with adapters (model-dependent: power pack with Euro-connector)	1
	Instructions for use, not shown	1

5.2 Assembling the device


Setting up the railing

- 1. Place the scale on a solid, level surface.
- 2. Release the latch for the railing.
- 3. Fold up the railing and hold it in this position.
- 4. Tighten up both locking screws until the railing is firm without any play.

Assembling the display housing

- 1. Take the display cable out of the fixing in the railing.
- 2. Connect the display cable in the railing to the cable of the display housing.

- 3. Push the cable carefully into the railing.
- 4. Place the display housing on the railing.
- 5. Ensure that the bores on the display housing and the bores in the railing are in line over one another.
- 6. Screw the display housing firmly to the railing.

Final steps

- 1. Ensure that all parts are correctly located.
- 2. Ensure that all screws are fully tightened.

5.3 Establishing power supply

Power is supplied to the scale with a power pack.

Connecting the power pack


WARNING!

Using the wrong power packs may cause bodily injury or damage to the device

Conventional power packs may deliver a higher voltage than is indicated on them. The scale may overheat, catch fire, melt or short-circuit.

- Use only original seca plug-in power packs with a 9 V or controlled 12 V output voltage.
- 1. If required, insert the power supply connector necessary for your power supply into the power pack.

- 2. Insert the connector of the power pack into the connecting socket of the scale.
- 3. Plug the power pack into a power supply socket.
- 4. Perform a function check as described in the section entitled "Function check" on page 89.

6. OPERATION


CAUTION!

Bodily injury

Perform a function check as described in the section entitled "Function check" on page 89 before using the device each time.

6.1 Making the scale operational

Moving the scale

- Disconnect the power pack connector from the power supply socket.
- 2. Disconnect the power pack connector from the connecting socket on the scale.
- 3. Undo the locking screws and carefully fold down the railing until it abuts completely.
- 4. Fold up the latch until it engages securely in the railing.
- 5. Fold away the access ramps.
- 6. Carefully align the device with the railing.
- Move the device in this position to the desired setup or storage location.


Aligning the scale


ATTENTION!

Incorrect measurement due to force shunts

If the scale and housing is in contact with something, e.g. a towel, weight will not be measured correctly.

- Set the scale up so that only its foot screws are in contact with the floor.
- 1. Place the scale on firm, level surface.
- 2. Undo the knurled wheels.


 Tighten the knurled wheels in the direction of the arrow.

The foot screws are secured against being adjusted.

Swiveling the display housing

The display housing of the device can be swiveled. As a result, it can be perfectly positioned for every application.

 Swivel the display housing so that you can operate and read it conveniently.

6.2 Weighing

CAUTION!

Patients can be injured if they fall

People with restricted mobility may fall if they try to step onto the scale or sit down on the scale.

 Support people with restricted mobility when they step onto the scale and when they sit down

Switching on the scale

Press the Start key.

All the elements of the display are shown briefly, then **SECA** appears in the display.

The scale is operational when **0.00** appears in the display.

Starting the weighing procedure

- 1. Ensure that the scale has no load.
- 2. Ask the patient to step onto the scale.
- 3. Ask the patient to keep still.
- 6350^{kg}

4. Read off the measurement result.

Determining extra weight (TARE)

hold

tare

hold

tare

Using the TARE function, you can prevent any extra weight (e.g. a towel or a support on the weighing platform) from affecting the weighing result.

ATTENTION!

Incorrect measurement due to force shunts

If an extra weight, e.g., a large towel, contacts the surface on which the scale is placed, weight will not be measured correctly.

- Make sure that extra weights are only placed on the scale's weighing surface.
- 1. Switch on the scale.
- 2. Place the extra weight on the scale.
- 3. Hold the (hold/tare) arrow key down until the message "NET" appears in the display.
- 4. Wait until the display stops flashing and 0.00 appears instead.
- 5. Weigh the patient as described in the section entitled "Starting the weighing procedure".
- 6. Read off the measurement result. The additional weight is automatically deducted.
- 7. To deactivate the TARE function, press the (hold/ tare) arrow key until the message "NET" is no longer displayed or turn off the scale.

NOTE


The maximum weight which can be displayed is reduced by the weight of the objects already placed on the scale.

Keeping the measurement result in the display (HOLD)

When you activate the HOLD function, the weight value continues to be displayed after the weight has been removed from the scale. As a result, you can attend to the patient before recording the weight.

- Ensure that the scale has no load.
- 2. Switch on the scale.
- 3. Weigh the patient as described in the section entitled "Starting the weighing procedure".
- 4. Press the (**hold/tare**) arrow key briefly.

The display flashes until a stable weight is measured. The weight is then continuously displayed. The \triangle symbol (non-calibratable function) and the message "HOLD" are displayed.

To deactivate the HOLD function, briefly press the (hold/tare) arrow key.

The \triangle symbol and the "HOLD" message are no longer displayed.

NOTE

If the auto hold function is activated, the weight value is automatically displayed permanently until the scale switches itself off or is switched off, see "Activating the auto hold (Ahold) function" on page 78.

Determining and rating body mass index (BMI)

Body mass index relates height and weight to one another and consequently allows more accurate information than ideal Broca weight, for example. A tolerance range is quoted which is considered ideal for health.

The device has three memories for height. You can enter and save the height of specific patients. Alternatively, you can save different start values and thus set a patient's actual height more quickly.

- 1. Ensure that the scale has no load.
- 2. Switch on the scale.
- Press the (bmi/menu) arrow key briefly.
 The message "BMI" appears.
 The last memory used is displayed (in this case, memory 2).
- 4. You can adopt the memory displayed or select a different memory using the arrow keys.
- 5. Confirm your setting with Enter (**send/print**). Arrows flash in the display.

The last height saved to the selected memory is displayed.


6. You can accept the height displayed or use the arrow keys to set a different height.


Confirm your setting with Enter (send/print).
 The height entered is saved and will be available again for the next BMI calculation.

NOTE

Make a note of the memory number so that you can call up this height for another BMI calculation.

 Weigh the patient as described in the section entitled Starting the weighing procedure.
 The patient's BMI is automatically calculated and displayed.

send

print

- 9. Read off the BMI and compare it with the categories quoted further down.
- 10. To deactivate the BMI function, briefly press the Enter key (**send/print**).

ВМІ	Rating		
Below 18.5	Patient is underweight. There could be anorexic tendencies. Weight should be increased to improve well-being and performance. If in doubt, seek the advice of a consultant.		
Between 18.5 and 24.9	The patient's weight is normal.		
Between 25 and 30 (incipient obesity)	The patient is slightly to moderately overweight. He or she should lose weight if already suffering from a disease such as diabetes, high blood pressure, gout or disorders of fat metabolism.		
Over 30	Weight needs to be lost urgently. The metabolism, the circulatory system and the bones are all under strain. A strict diet, lots of physical activity and behavioral training are all recommended. If in doubt, seek the advice of a consultant.		

Sending measurement results to wireless receivers

If the scale is integrated in a **seca 360° wireless** network, you can send the measurement results to devices ready to receive them (e.g. wireless printer, PC with USB wireless adapter) at the touch of a button.

- Press the Enter key (send/print).
 - Press key briefly: send measurement results to all devices ready to receive them
 - Press and hold key: print out measurement result on wireless printer

Calculating and printing out BMI automatically

If you use this scale together with a wireless printer and a length measuring rod from the **seca 360° wireless** system, you can have BMI calculated and printed out automatically.

NOTE

The prerequisite for this function is that the devices are registered together in one wireless group (see "The seca 360° wireless network" on page 82).

- 1. Perform the weighing operation.
- Briefly press the Enter key (send/print) of the scale.

The measured value is sent to the wireless printer but not printed.

- 3. Perform the height measurement.
- 4. Press and hold the Enter key (**send/print**) of the length measuring rod.

The measured value is sent to the wireless printer. BMI is calculated.

Height, weight and BMI are printed out.


Switching off the scale

Press the Start key.

6.3 Other functions (menu)

Other functions are available to you in the menu for the scale. This enables you to configure the scale to suit your conditions of use perfectly.

*The description of menu item "rF" can be found in the section entitled "Operating the scale in a wireless group (menu)" on page 84.


Navigating in the menu

- 1. Switch on the scale.
- Keep arrow key (bmi/menu) pressed until the menu is called up.

The last menu item selected (in this case, auto hold "Ahold") appears in the display.

Press one of the arrow keys until the desired menu item appears in the display (in this case, filtering "Fil").

- Confirm your selection with Enter (send/print).
 The current setting for the menu item or a submenu are displayed (in this case, level "0").
- 5. To amend the setting or call up a different submenu, press one of the arrow keys until the desired setting (in this case, level "2") is displayed.
- 6. Confirm the setting with Enter (**send/print**). You exit the menu automatically.
- 7. To make more settings, call up the menu again and proceed as described.

NOTE

If no key is pressed for about 24 seconds, you exit the menu automatically.

Deleting saved values automatically (ACIr)

₽₽₽

Πu

To prevent current measurement results remaining in the memory of the device and leading to an incorrect BMI calculation, you can set the scale so that measurement results are automatically deleted after 5 minutes.

NOTE

This function is activated at the factory on some models. You can deactivate the function if desired.

- 1. Select the item "ACIr" from the menu.
- 2. Confirm your selection.
- 3. Select the desired setting:
 - or
 - off
- 4. Confirm your selection.

You exit the menu automatically.

Permanently saving the additional weight (Pt)

Using the Pre-tare function, you can permanently save an additional weight and automatically deduct it from a measurement result. For example, you can save an all-inclusive weight for shoes and clothing and then always deduct this from the measurement result if a patient is weighed fully-clothed.

The device has three memories for weights. You can save various weights and call them up individually depending on the initial situation, so that they are automatically deducted from the measurement result.

PE , '.

Select the item "Pt" from the menu.
 The message "Pt" appears.

2. Confirm your selection.

The last memory used is displayed.

- You can adopt the memory displayed or select a different memory.
- 4. Confirm your selection.

Arrows flash in the display.


The additional weight saved to the selected memory is displayed.

5. You can adopt the saved value or amend it using the arrow keys.


NOTE

When you enter the value "0", the function is switched off. The message "Pt" is no longer displayed.

- 6. Confirm your selection.
- 7. Ask the patient to position him or herself on the scale.

The saved additional weight has been deducted automatically.

9. Confirm your selection.

The function is deactivated.

You exit the menu automatically.

NOTE


When you switch the scale off, the function is switched off. The message "Pt" is no longer displayed when you switch the scale on again.


Activating the auto hold (Ahold) function

When you activate the auto hold function, the measurement result for every weighing operation continues to be displayed after the weight has been removed from the scale. It is then no longer necessary to activate the Hold function manually for each individual weighing operation.

NOTE

This function is activated at the factory on some models. You can deactivate the function if desired.

- 3. Select a filtering level.
 - 0: no filtering
 - 1: moderate filtering
 - 2: high filtering
- 4. Confirm your selection.

You exit the menu automatically.

Switching unit of weight (Unit)


On non-calibrated scales, you can select the unit (Unit) in which you want weight to be displayed.

CAUTION!

Patient hazard

In order to avoid misinterpretations, test results for medical use must be displayed and used in SI units only. Some devices offer the ability to display test results in other units. This is only an additional function.

- Use the results exclusively in SI units.
- The use of measurement results in non-SI units is the sole responsibility of the user.
- 1. Select the item "Unit" from the menu.
- 2. Confirm your selection.

Un 1E

The current setting is displayed.

Select the unit in which you want weight to be displayed:

- kilograms (kg)
- pounds (lbs)
- stones (sts)
- 3. Confirm your selection.

You exit the menu automatically.

Restoring factory settings (RESET)

You can restore the factory settings for the functions below.


Function	Factory setting
Auto hold (Ahold")	Model-dependent
Beep (Press)	Off
Beep (Hold)	On
Filtering (Fil)	0
Autoclear (Aclear)	Model-dependent
Pre-tare (Pt)	0 kg
Height for Body Mass Index (BMI)	170 cm

Function	Factory setting
Unit of weight	kg
Wireless module (SYS)	Off
Autosend (ASend)	Off
Autoprint (APrt)	Off

NOTE

The wireless module is switched off when factory settings are restored. Information about existing wireless groups is retained. Wireless groups do not have to be set up again.

- 1. Select the item "Reset" from the menu.
- 2. Confirm your selection.
 You exit the menu automatically.
- Switch off the scale.
 Factory settings are restored and are available when the scale is switched back on.

7. THE SECA 360° WIRELESS NETWORK

7.1 Introduction

This device is equipped with a wireless module. The wireless module allows wireless transmission of measurement results for evaluation and documentation. Data can be transmitted to the following devices:

- seca wireless printer
- PC with USB wireless adapter

seca wireless groups


The **seca 360° wireless** network operates with wireless groups. A wireless group is a virtual group of transmitters and receivers. If several transmitters and receivers of the same type are going to be operated, up to 3 wireless groups (0, 1, 2) can be set up with this device.

Setting up several wireless groups ensures the reliable transmission of correctly-addressed measured values if you wish to use several examination rooms each with comparable equipment.


The maximum distance between transmitters and receivers is approx. 10 meters. Certain local circumstances such as the thickness and type of walls, may reduce the range.

The following combination of devices is possible for each wireless group:

- 1 baby scale
- 1 personal scale
- 1 length measuring rod
- 1 seca wireless printer
- 1 PC with seca USB wireless adapter

Channels

Within a wireless group, devices communicate with each other on three channels (C1, C2, C3). This guarantees that data are transmitted reliably and without interference.

If you set up a wireless group with this scale, the device suggests three channels which will ensure optimal data transmission. We recommend adopting the channel numbers suggested.

You can also manually select the channel numbers (0 to 99), for example if you wish to set up several wireless groups.

To ensure interference-free data transmission, the channels must be spaced sufficiently far apart. We recommend a spacing of min. 30 between channel numbers. Each channel number may only be used for one channel.

Example configuration for channel numbers when setting up 3 wireless groups within a medical practice:

- wireless group 0: C1=_0, C2= 30, C3=60
- wireless group 1: C1=10, C2=40, C3=70
- wireless group 2: C1=20, C2=50, C3=80


Device detection

If you set up a wireless group with the scale, it searches for other active devices from the **seca 360°** wireless system. The devices detected are shown in the display of the scale in the form of modules (e.g. MO 3). The numbers have the following meaning:

- 1: personal scale
- 2: length measuring rod
- 3: wireless printer
- 4: PC with seca USB wireless adapter
- 7: baby scale
- 5, 6 and 8-12: reserved for system expansion

Operating the scale in a wireless group (menu)

All the functions you need to operate the device in a seca wireless group can be found in the "rF" submenu. Information for navigating in the menu can be found on Page 76.

Activating the wireless module (SYS)


The device is supplied with the wireless module deactivated. You have to activate it before you can set up a wireless group.

- 1. Switch on the device.
- 2. Select the "SYS" menu item from the "rF" submenu.
- 3. Confirm your selection.
- 4. Select the "on" setting.
- 5. Confirm your selection. You exit the menu automatically.

group (Lrn)

To set up a wireless group, follow the instructions below.

- 1. Switch on the device.
- 2. Call up the menu.
- 3. Select the item "rF" from the menu.
- 4. Confirm your selection.
- 5. Select the "Irn" (learn) menu item from the "rF" submenu.
- 6. Confirm your selection.

19 1		If wireless group "0" already exists, use the arrow keys to select a different ID (in this case: wireless group 1 "ID 1").
	7.	Confirm your selection of wireless group.
		The scale suggests a channel number for Channel 1 (in this case, C1 "0").
		You can adopt the proposed channel number or use the arrow keys to set a different channel number.
	8.	Confirm your selection for Channel 1.
[230]		The device suggests a channel number for Channel 2 (in this case, C2 "30").
		You can adopt the proposed channel number or use the arrow keys to set a different channel number.
		NOTE Double-digit channel numbers are displayed without a space. The display "C230" therefore

The wireless group currently set (in this case:

wireless group 0 "ID 0") is displayed.

5£0P

r 360

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10. Confirm your selection for Channel 3.

9. Confirm your selection for Channel 2.

3 (in this case, C3 "60").

The message **StOP** appears in the display. The device waits for signals from other wireless-capable devices in range.

means: Channel "2", channel number "30"

The device suggests a channel number for Channel

You can adopt the proposed channel number or use the arrow keys to set a different channel

NOTE

number.

On some devices, a special switch-on procedure has to be followed if they are going to be integrated in a wireless group. Follow the instructions for use for the device in question.

11. Switch on the device, e.g. a wireless printer, you want to integrate in the wireless group.

When the wireless printer is detected, you will hear a beep.

NOTE

As soon as you have integrated a wireless printer in the wireless group, you then have to select a print option (Menu\rf\APrt) and set the time (Menu\rf\time).

- 12. Repeat step 11. for all the devices you want to integrate in this wireless group.
- 13. Press the Enter key to stop the search procedure.

14. Press one of the arrow keys to have the devices which have been detected displayed (in this case, Mo 3 for a wireless printer).

If you have integrated several devices in one wireless group, press the arrow keys several times to ensure that all devices have been detected by the scale.

15. Exit the menu using the Enter key or wait until you exit the menu automatically.

Activating automatic send (ASend)

You can configure the device so that the measurement results are automatically sent to all receivers which are ready to receive and are registered to the same wireless group (e. g.: wireless printer, PC with USB wireless adapter).

NOTE

If you are using a wireless printer, ensure that "off" is not set as the print option (see "Selecting print option (APrt)" on page 86).

- 1. Switch on the device.
- ASEnd
- **□**n
- 2. Select the "ASend" menu item from the "rf" submenu and confirm your selection.
- 3. Select the "on" setting and confirm your selection. You exit the menu automatically.

Selecting print option (APrt)

You can configure the device so that the measurement results are automatically printed out on a wireless printer registered to the wireless group.

NOTE

This function is only accessible if a seca wireless printer has been integrated in the wireless group via the "learn" function.

1. Switch on the device.

 $\Pi \Pi$

- 2. Select the "APrt' menu item from the "rf" submenu and confirm your selection.
- Select the setting relevant for your device combination:
 - HI: measurement results from length measuring devices
 - MA: measurement results from scales
 - HI_MA: measurement results from length measuring devices and scales
 - off: no automatic printing, printing only possible by pressing and holding the Enter key during the weighing procedure.
- 4. Confirm your selection.

You exit the menu automatically.

Set the time (Time)

You can configure the system so that the wireless printer automatically adds a date and time to your measurement results. To do this, you have to set date and time once on this device and transmit it to the internal clock of the wireless printer.

NOTE

This function is only accessible if a seca wireless printer has been integrated in the wireless group via the "learn" function.

- Switch on the device.
- Select the "Time" menu item from the "rF" submenu.
- Confirm your selection.
 The current setting for "Year" is displayed.
- 4. Set the correct year.
- 5. Confirm your selection.
- Repeat steps 3. and 4. accordingly for "month" (Mon), "day" (dAy), hour (hour) and minute (Min).
- 7. Confirm each of your selections.

After confirming the setting for minute, you automatically exit the menu.

The settings are automatically transmitted to the wireless printer.

The wireless printer automatically adds a date and time to every printout.

NOTE

Follow the instructions for use for the wireless printer for further information about operating it.

8. HYGIENE TREATMENT

WARNING! Electric shock

The device is not de-energized when the on/off key is pressed and the display goes out. Use of fluids on the device may cause electric shock.

- Before each hygiene treatment, remove the mains cable to de-energize the device.
- Ensure that no fluids penetrate the device.

CAUTION!

Damage to device

Unsuitable cleaning agents may damage the sensitive surfaces of the device.

- Do not use aggressive or abrasive cleaning agents.
- Do not use white spirit or petroleum spirit.

8.1 Cleaning

 Use a soft cloth dampened with mild soapsuds to clean the surfaces of the device.

8.2 Disinfecting

- 1. Follow the instructions on the disinfectant.
- Disinfect the device at regular intervals with a commercially available disinfectant, e.g. with 70% isopropyl alcohol.

Please take note of the following terms:

- Prior to every measurement with direct skin contact:
 - Weighing platform
- After every measurement with direct skin contact:
 - Weighing platform
- If required:
 - Display
 - Membrane keyboard

8.3 Sterilizing

The sterilisation of this device is not permitted.

9. FUNCTION CHECK

♦ Perform a function check before each use.

A complete function check includes:

- a visual inspection for mechanical damage
- a visual inspection and function check of the display
- function check of all the controls shown in the section entitled "Overview".
- a function check of the optional accessories

If you find faults or deviations during the function check, as a first step try to remedy the fault with the aid of the section entitled "Troubleshooting" from page 89.

CAUTION!

Bodily injury

If you find faults or deviations during the function check which you are unable to remedy with the aid of the section entitled "Troubleshooting" from page 89, you may not use the device.

 Have the device repaired by seca service or an authorized specialist dealer.

10.TROUBLESHOOTING

Problem	Cause / Solution	
No weight displayed with a load on the scale.	The scale has no power supply. - Check whether the scale is switched on.	
0.00 does not appear before weighing.	A load was already on the scale before it was switched on. - Remove the load from the scale. - Switch the scale off and back on again.	
One segment is illumi- nated either continu- ously or not at all.	The corresponding point has a fault Inform the service department.	
The display "StOP" appears.	Maximum load has been exceeded Remove the load from the scale.	
The display "tEMP" appears.	The ambient temperature of the scale is too high or too low. - Set up the scale in an ambient temperature between +10 °C and +40 °C. - Wait for around 15 minutes until the device has adapted to the ambient temperature.	

Problem	Cause / Solution		
After switching on for the first time, measurement results are transmitted and two beeps are heard.	 The device was unable to send measurement results to the wireless receiver (seca wireless printer or PC with seca USB wireless adapter). Ensure that the scale is integrated in the wireless network. Ensure that the receiver is switched on. Reception is being interfered with by HF equipment (e.g. cell phones) in the vicinity. Keep HF equipment at least 1 meter away from transmitters and receivers in the seca wireless network. The actual transmission output of HF equipment may require minimum distances of more than 1 meter. For details, go to www.seca.com. 		
	NOTE If this fault is not eliminated, there will be no fresh acoustic warning when further attempts are made to send data.		
Only the "SYS" item is visible in the rf menu.	The wireless module is deactivated. Activate wireless module (see "Activating the wireless module (SYS)" on page 84).		
Only the "SYS" and "Irn" items are visible in the rf menu.	 The wireless module is activated and no wireless group is set up. Set up wireless group (see "Setting up a wireless group (Lrn)" on page 84). 		
The items "APrt" and "Time" are not visible in the rf menu.	 No wireless printer registered to the wireless group. Register the wireless printer to the wireless group via the "Irn" menu item (see "Setting up a wireless group (Lrn)" on page 84). 		
After calling up the menu, the "rf" item is not displayed.	The scale's wireless module is defective. Inform seca service.		
The display "Er:X:11" appears.	The scale has too high a load or too high a load in one corner. - Take the load off the scale or distribute the weight more evenly. - Restart the scale.		
The display "Er:X:12" appears.	The scale has been switched on with too high a load. - Remove the load from the scale. - Restart the scale.		

Problem	Cause / Solution		
The display "Er:X:16" appears.	The scale was caused to oscillate, the zero point could not be determined Restart the scale.		
The Enter key (send/ print) is pressed and the display "Er:X:71" appears.	No data transmission possible, wireless module deactivated. - Activate wireless module (see "Activating the wireless module (SYS)" on page 84).		
The Enter key (send/ print) is pressed and the display "Er:X:72" appears.	No data transmission possible, no wireless group set up Set up wireless group (see "Setting up a wireless group (Lrn)" on page 84).		

11.SERVICING

On leaving the factory, your seca scale has an accuracy of $\pm 0.15\%$ or better. To preserve this level of accuracy, the product must be set up with care and serviced regularly. Depending on how frequently the scale is used, we recommended servicing at intervals of 3 to 5 years.

ATTENTION!

Incorrect measurements as a result of poor servicing

- Please only have servicing and repairs performed by an authorized service partner.
- You can find service partners in your area at www.seca.com or by sending an e-mail to service@seca.com.

12.TECHNICAL DATA

Technical data seca 676			
Dimensions of scale when folded up Depth Width Height	965 mm 1089 mm 151 mm		
Dimensions of scale when set up Depth Width Height	1150 mm 920 mm 1120 mm		
Dimensions of weighing platform • Depth • Width • Height	970 mm 800 mm 54 mm		
Net weight	approx. 34.7 kg		
Temperature range	+10 °C to +40 °C -10 °C to +65 °C -10 °C to +65 °C		
Air pressure Operation Storage Transport	700-1060 hPa 700-1060 hPa 700-1060 hPa		
Humidity Operation Storage Transport	30% - 80% no condensation 0% - 95% no condensation 0% - 95% no condensation		
Height of digits	25 mm		
Power supply	Power pack		
Power consumption • with wireless module deactivated • with wireless module activated	approx. 20 mA approx. 37 mA		
Medical device in accordance with Directive 93/42/EEC	Class I with measuring function		
EN 60601-1: • insulated device, protection class II			
electrical medical device, type B	<u>★</u>		
Maximum load	360 kg / 800 lbs / 57 sts		
Minimum load	1.0 kg		
Increments	50 g / 0.1 lbs		

Technical data seca 676			
Tare range	360 kg / 800 lbs / 57 sts		
Accuracy • 0 to 67 kg • 67 kg to maximum load • 0 to 75 lbs • 75 lbs to maximum load • 0 to 10.55 sts • 5.5 sts to maximum load	± 100 g ±0.15% ± 0.2 lbs ±0.15% ± 0.2 lbs ±0.15%		
Wireless transmission Frequency band Transmission power Standards applied	2.433 GHz - 2.480 GHz < 10 mW EN 300 328 EN 301489-1: EN 301489-17:		

13.OPTIONAL ACCESSORIES

seca 360° wireless devices	Article number		
Length measuring rods			
• seca 274	Country-specific versions		
• seca 264	Country-specific versions		
Access ramp	470-00-00-009		
Wireless printers			
seca 360° Wireless Printer 465	Country-specific versions		
seca 360° Wireless Printer Advanced 466	Country-specific versions		
PC software			
seca analytics 115	Application-specific license		
	packages		
seca 360° Wireless USB Adapter 456	456-00-00-009		

14.SPARE PARTS

Spare parts	Article number
Power pack with Euro-connector: 230 V~ / 50 Hz/12 V= / 130 mA	68-32-10-252
Switchmode power pack with adapters: 100 - 240 V~ / 50 -60 Hz / 12 V= / 0.5 A	68-32-10-265

15. DISPOSAL

Do not dispose of the device with household waste. The device must be disposed of properly as electronic waste. Comply with the national provisions applicable in your country. For further information contact our service department at:

service@seca.com

16.WARRANTY

We offer a two-year warranty from the date of delivery for defects attributable to faulty material or poor work-manship. This excludes all moveable parts such as (rechargeable) batteries, cables, power supply units, etc. Defects which are covered by the warranty shall be rectified free of charge for customers on production of the sales receipt. No further claims can be accepted. The costs of shipment in both directions shall be borne by the customer where the device is not located at the customer's premises. In the event of any damage during shipment warranty claims can only be asserted where the complete original packaging was used for shipment and the scales were secured inside in the same manner as in the original packaging. You should therefore keep all packaging.

The warranty shall become null and void where the device is opened by persons not expressly authorised to do so by seca.

We ask customers based abroad to contact their local sales agent direct in the case of warranty claims.

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FCC ID: X6T172A01

IC: 8898A-172A01

This device complies with Part 15 of the FCC Rules and with RSS-210 of Industry Canada. Operation is subject to the following two conditions. (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE:

This device complies with Part 15 of the FCC Rules and with RSS-210 of Industry Canada. Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

NOTE:

Changes or modifications made to this equipment not expressly approved by seca may void the FCC authorization to operate this equipment.

NOTE:

Radiofrequency radiation exposure Information: This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance of 1 m between the radiator and your body. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

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The scales comply with the applicable requirements of the following directives:

93/42/EEC governing medical devices.

Furthermore the following directives are applicable among others:

FN 300 328, FN 301 489-1 and -17 governing electromagnetic compatibility and radio spectrum

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