
seca 763

D	<i>Bedienungsanleitung und Garantieerklärung</i>	4
GB	<i>Instruction manual and guarantee</i>	18
F	<i>Mode d'emploi et garantie</i>	31
I	<i>Manuale di istruzioni e garanzia</i>	45
E	<i>Manual de instrucciones y garantía</i>	59
P	<i>Instruções de utilização e declaração de garantia</i>	73

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
1. Signed and sealed

With seca products, you are not only purchasing technology which has been perfected over 100 years, but also quality with official legal certification confirmed by institutions. seca products comply with European directives, standards and national laws. When you buy seca, you buy the future.

Products bearing this symbol meet the applicable regulatory requirements of the European Community, especially the following:

- Directive 93/42/EEC for medical devices

seca's professionalism is also recognized by official testing agencies. TÜV Süd Product Service, the appointed office for medical devices, confirms with this certificate that as a medical devices manufacturer, seca consistently complies with the strict legal requirements. seca's quality assurance system includes the areas of design, development, production, sales and service of medical scales and length measuring systems as well as software and measuring systems for assessing state of health and nutrition.

seca helps the environment. Saving natural resources is very important to us. We therefore make every effort to save on packaging materials wherever it makes sense and whatever is left over can be conveniently disposed of on site via the recycling system.

2. Congratulations!

Your purchase of this integrated measurement station **seca 763** means you have a high-precision and sturdy electronic body height and body weight measuring instrument. seca has been putting its experience at the service of health for over 150 years now, and as market leader in many countries of the world, is always setting new standards with its innovative developments for weighing and measuring.

The integrated **seca 763** measurement station is primarily used in accordance with national regulations in hospitals, doctor's practices and health care centers.

The integration of electronic scales and electronic height measurement in one device means the **seca 763** is particularly

suitable where weighing and measuring must be carried out where little space is available. The data measured can be converted at the push of a button into a Body Mass Index. The telescoping design used for the height measurement unit as well as the retractable head slide allows the device to be easily moved on the built-in castors.

The sturdy design and intelligent power management system ensure a long period of maintenance-free service so that the device can be used for years to come. Operation is easy with two large displays allowing weight and height to be read out simultaneously.

3. Safety

Before using the device, please take some time to read the following safety instructions.

- Follow the safety instructions in the instruction manual.
- Keep the operating manual and the declaration of conformity contained therein in a safe place.
- Never drop the device or expose it to strong shocks.
- Have repairs carried out exclusively by authorised persons.
- After moving the device, place it carefully back in the position for use.
- After using the height measurement unit fold the head slide down. This avoids the hazard of injury.
- Be sure to use only the seca AC adaptor supplied. Before plugging the AC adaptor into a socket, check that the mains voltage indicated on the AC adaptor matches local mains voltage.
- When using the device with an AC adaptor, be sure that the power supply cable does not pose any danger of tripping or entanglement.
- Be sure to have the device maintained regularly (see "Maintenance").


4. Before you start...

Unpacking

Remove the packaging and place the based on a solid, even surface.


The equipment supplied includes:

- Base of scales
- Column with integrated height measurement and display head
- Protective cover for column
- SW 3 Allen spanner
- SW 5 Allen spanner
- 2 screws
- seca AC adaptor
- Operating instructions with a declaration of conformity

Mounting the column

- Place the protective cover onto the column.
- Place the column in to the notch so that the folding head slide is pointing down toward the platform.
- Screw the column from beneath onto the weighing platform.
- If needed, carefully place the scales on the side to tighten the screws.
- As shown in the illustration, connect the power supply cable on the bottom of the weighing platform.
- Insert the plug of the display cable in the jack on the bottom of the weighing platform.
- Fasten the cable in the cable guide.
- Now put the device back in an upright position.
- Slide the column cover all the way down until it locks into place on the weighing platform.


Power supply

The power supply must come from the supplied AC adaptor.

The connector for the AC adaptor is located under the rod.


Please only use the power supply unit supplied with the product or power supply units from the seca accessories range. Before plugging into the socket, check whether the mains voltage details on the power supply unit match those of the local voltage supply.

Please use only the mains unit supplied or mains units from the seca range of accessories. Before plugging the unit into the socket, check whether the mains voltage marked on the mains unit matches local mains voltage.

5. Setting up and aligning the scales


- Place the scales on a level surface.
- Align the scales by removing the foot screws. The scales must stand securely on the floor on all four feet.

6. Operation

Controls and displays

The controls and display are in the display part.

START

Switching the device on and off

kg lbs sts

Press button briefly → switch between kg and lbs

hold down → switch between kg and sts

In the function menu → Raise the value

HOLD TARE

Press button briefly → Active/deactivate hold function

hold down → Set the display to zero, activate/deactivate TARE

In the function menu → Lower the value

FUNC SET

Activate BMI function, confirm

cm / inch

Toggling the unit of height measurement between cm and inch

ZERO

Setting the height measurement to zero for relative measurements

DATA VALID

Pushing button to mark the valid height measurement value

Correct weighing

- Press the Start key with no load on the scale.
SECA, **88.888**, and **0.0** appear consecutively in the display. The scale is then automatically set to zero and ready to use.
- If necessary, change over the weight display (see "Switching the weight display").
- Stand on the scales and remain motionless.
- Read the weighing result off the digital display.
- If a weight of more than 250 kg or 550 lbs is placed on the scales, then **STOP** appears in the display.
- To switch off the scales, press the Start button again.

Switching the weight display

The weight display can be switched between kilogrammes (kg) and pounds (lbs) or between kilogrammes (kg) and stones (sts).

- Switch on the scale using the green **START** key. Weight is initially displayed in the last setting selected.
- To switch the weight display between kilogrammes and pounds, press the **kg/lbs/sts** key.
- The **kg** or **lbs** display lights up.

- To switch weight display between pounds and stones, hold down the **kg/lbs/sts** key for approx. 2 seconds.
- The **lbs** or **sts** display lights up.

Switching off the device

- To switch off the scales, press the Start button again.

Taring weight display (TARE)

The TARE function is useful if an additional weight should not be considered while weighing.

Proceed as follows:

- Press the Start key with no load on the scale.
- First, place the additional weight onto the scales and activate the TARE function by pressing the **HOLD/TARE** button.

The display is reset to zero. **0.0** appears in the display. The note "NET" appears.

- Now place the patient onto the scales with the previously weighed additional weight.

The scales show the weight of the patient without the additional weight. You can now perform as many weighing operations as you like.

- Deactivate the TARE function by pressing the button again; this puts the scales back into the normal weighing mode, the tare value (the additional weight value) is deleted.
- The tare value is likewise erased when the scales are switched off.

Store weight value (HOLD)

You can also continue displaying the weight determined once the load has been removed.

This makes it possible to look after a patient first, before having to note the weight:

- Press the **HOLD/TARE** button while the scales have a weight on it.
- "HOLD" appears in the digital display. The weight remains stored in the display until the Hold function is automatically switched off.
- The triangle with the exclamation mark indicates that this frozen display is an additional facility which is not subject to calibration.
- The HOLD function can be switched off by pressing the **HOLD/TARE** key again. **0.0** appears in the display if the scales have no load, otherwise the current weight is shown. The "HOLD" display vanishes.

Note: The HOLD and TARE functions can also be used simultaneously.

Height measurements

- Press the **START** button on the display housing to switch on the device.

First ----- appears in the height display. The unit selected before switching off the device is active and is marked by an arrow.


- If required, switch the unit of height measurement (see "Toggling unit of height measurement").
- Raise the head slide up and push the head slide up far enough for the person whose height is being measured to easily stand under it.
- Ask the person being measured to place his back against the height rod. Make sure the person's posture is upright and head is straight.
- Now push the head slide down until it rests on the person's head.
- Read off the height on the display.

Toggling unit of height measurement

The display can toggle between "cm" and "ins" (inches).

- Switch the device on. The height read-out is displayed in the most recent unit selected.
- To switch the unit displayed press the **cm/ins** button.
- An arrow lights up next to the **cm** or **ins** display.

Relative measurements

The height rod can be set to zero at any position. This function allows extremities to be measured.

- Press the **ZERO** button in order to set the height rod to zero at its current position.

ZERO SET is highlighted in the display.

All lengths are now measured starting from this zero set point. Any measurement below this zero point will cause results to be displayed with a minus sign.

- Press the **ZERO** button again if you want to start measuring from absolute zero again.

Body Mass Index (BMI)

The Body Mass Index compares body size with body weight while also allowing more precise information to be determined for ideal body weight according to Broca. A tolerance range is given which is considered ideal for health.

By using the integrated measurement station **seca 763** you can have the Body Mass Index calculated directly from the size measured and the measured weight. Please follow these steps.

- Press the Start key with no load on the scale.

SECA, **8.8.8.8.8** and **0.0** appear consecutively in the display. The scale is then automatically set to zero and ready to use.

- Now measure height as described. After you have determined that the patient is standing up straight and the head slide is correctly positioned, then please push the **DATA VALID** button.

This marks the height measurement as a valid value for the Body Mass Index calculation.

- Push the head slide back up.

- Now read the weight value in the weight display.

- Press the FUNC key. The measured height which is marked as valid now appears in the display. Push the **kg/lbs/sts** buttons (raise the value) or **HOLD/TARE** (lower the value) if you want to correct the height, for example to compensate for shoe heel height.

- Push the **FUNC** button again. The BMI function is now activated.
- Now the patient should try to stand as still as possible.
- Read the BMI on the LED and compare it with the categories under "Evaluating Body Mass Index".
- To return to normal weighing mode, press the **FUNC** key.
- To switch off the scales, press the Start button again.

Note:

The Body Mass Index is correctly calculated even when you have not selected metric display units.

Evaluating Body Mass Index for adults:

Compare the value determined with the categories quoted below which correspond to those used by WHO, 2000.

- **BMI below 18.5:**

Warning, this patient weighs too little. There could be a tendency to anorexia. An increase in weight is recommended to improve well-being and performance. If in doubt, consult a specialist.

- **BMI between 18.5 and 24.9:**

The patient is of a normal weight. He can stay as he is.

- **BMI between 25 and 30 (pre-obese):**

The patient is slightly to moderately overweight. The patient should lose weight if a previous illness exists (such as diabetes, high blood pressure, gout, lipid disorders).

- **BMI over 30:**

Weight reduction is essential. Metabolism, cardiovascular system and bones are all under strain. A consistent diet, plenty of exercise and behavioural training are all recommended. If in doubt, consult a specialist.

7. Care and cleaning

When needed, clean the device with a household cleaning agent or other common disinfectant. Follow the manufacturer's instructions.

Under no circumstances use abrasive or acid cleaners, white spirit, benzene or the like for cleaning. Such substances can damage the high-quality surfaces.

8. What to do if ...

... no weight or height readout appears?

- Is the device switched on?
- Is the plug between the power supply connector in the rod and base properly connected?
- Is the plug contact of the AC adaptor properly connected on the device rod?
- Is the plug contact for connecting the display and the base properly connected?

... 0.0 does not appear before weighing?

- Press the start button again – there must not be any load on the scales – and only its feet should be in contact with the floor.

... one segment is illuminated constantly or not at all?

- The relevant segment has a fault. Inform Servicing.

... STOP appears in the weight display?

- Maximum load has been exceeded.

... the weight display is blinking?

- If you have just activated the Hold or Tare function, please wait until the scales have detected a steady weight value. The display will then stop flashing.
- If there has been no activation of a function beforehand, remove the load from the scales and wait until 0.0 is displayed, then weigh again.

... the display EIP appears?

- The ambient temperature of the device is too high or too low. Place the device where the ambient temperature is between +10 °C and +40 °C. Wait about 15 minutes until the device has adjusted to the ambient temperature and weigh again.

... the display E and a number appear?


- Use the Start button to switch off the device and start again. After this the device will function normally.
- If this is not the case, disconnect the power supply by briefly unplugging the AC adaptor. If this measure is equally unsuccessful, inform the Service department.

... ----- appears in the height display?

- Slide the head slide up or down quite a ways. This will make the head slide pass a mark which can be evaluated and the system can now calculate the position correctly again.

... the head stop slides down and does not stay in the required position?

- Place the scale carefully on its side so that the base of the scale is accessible. Insert the Allen key supplied in the rectangular hole on the end of the measuring rod as shown in the illustration and tighten up the screw for the measuring rod guide clockwise.
The measuring rod with the head stop should remain freely mobile, but not be able to slide down on its own.
Then stand the scale upright again.

9. Servicing

Your seca scale leaves the factory with an accuracy greater than $\pm 0,15\%$. To ensure that this accuracy is maintained, the product must be carefully set up and regularly maintained. Depending on the frequency of use, we recommend servicing every 3 to 5 years. Please ensure that a qualified service agent is used. If in doubt, contact seca Customer Service for further information.

10. Technical data

Device dimensions

Width:	360 mm / 14.2 ins
Height:	1320 +10 mm / 52 +0.4 ins
Depth:	686 mm / 27 ins

Weight approx 20.5 kg

Temperature range +10 to +40 °C

Insulated equipment, Protection Class II:

medical electrical equipment type B

Colour white/black

Maximum load 250 kg / 550 lbs / 39 sts

Graduations 50 g / 0.1 lbs

Precision better than $\pm 0.15\% / \pm 100\text{g}$

Height measurement

Precision Tolerance of $\pm 2\text{ mm}$

Reproducibility $\pm 1\text{ mm}$

Measurement range 110 - 200 cm / 43 - 79 ins

11. Spare parts

RS232 adapter kit with print port

Order no. 460-00-00-004

Model-specific power unit

Designation:

Plug-type AC adaptor 230 V~/50 Hz/12 V=/130 mA

Order no. 68-32-10-252

Designation:

Switchmode power adapter,
100-240V~ / 50-60Hz / 12V= / 0.5A

Order no. 68-32-10-265

12. Disposal

Disposing of the device

Do not dispose of the device in domestic waste. The device must be disposed of properly as electronic waste.

Follow the national regulations which apply in your case. For further information, contact our service department at:

service@seca.com

13. Warranty

A two-year warranty from date of delivery applies to defects attributable to poor materials or workmanship. All moveable parts - batteries, cables, AC adaptors, rechargeable batteries etc. - are excluded. Defects which come under warranty will be made good for the customer at no charge on production of the receipt. No further claims can be entertained. The costs of transport in both directions will be borne by the customer should the equipment be located anywhere other than the customer's premises. Warranty claims for

damage during transport can only be enforced if the complete original packaging was used in transit and the goods contained are secured and fixed into position in their original condition. All the packaging should therefore be retained.

A claim under warranty will not be honoured if the equipment is opened by persons not expressly authorised by seca to do so.

We would ask our customers abroad to contact their local sales agent in the event of a warranty matter.

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