seca analytics **115** Stručný návod k instalaci



Stručný návod 🤇 seca analytics 115





Stručný návod k instalaci Seřízení jednotlivého pracoviště

Tento stručný návod Vám ukáže, jak můžete rychle a snadno uvést do provozu Váš software seca analytics 115 společně s měřicími přístroji seca 360° wireless.

Potřebné jsou pouze 4 kroky:

1. Instalace softwaru	Str. 3-5
2. Aktivace licence jednotlivého pracoviště	Str. 6-8
3. Změna jazyka dialogu	Str. 9-10
4. Konfigurace bezdrátového provozu	Str. 11-14

Pokyny týkající se provádění měření, vyhodnocení výsledků vyšetření, provozu sítě, správy uživatelů a všech dalších funkcí softwaru viz podrobná uživatelská dokumentace skládající se z příručky administrátora a návodu k použití na DVD. Tento stručný návod Vám ukáže cestu také k těmto dokumentům. **Při instalaci potřebujete na Vašem počítači administrátorská práva.**

1. Instalace softwaru

Krok č. 1

Vložte DVD. **Pokyn:** Pokud se DVD nespustí automaticky, otevřete ve Windows Explorer složku DVD a dvakrát klikněte na "Setup".

😾 seca analytics 115 Setup	1
Seca analytics 115 Welcome Setup will install seca analytics 115 on your computer Click install to continue, options to set the install directory or Close to exit. Select option: analytics 115: complete Projects 115: server	
seca analytics 115 <u>license terms</u> . Version 1.4.0.0 Lagree to the license terms and conditions Qptions install Close	
↓	
🚽 seca analytics 115 Setup	1
seca analytics 115	
Setup will install seca analytics 115 on your computer. Click install to continue, options to set the install directory or Close to exit.	
Select option: analytics 115: complete analytics 115: client analytics 115: server	
seca analytics 115 <u>license terms</u> .	
Version 1.4.0.0	
Qptions Provide United to the Contract of the	



Postupovat podle pokynů na obrazovce





Po restartu sys	stému se zobra	zí automatick	ky:				
	Jednotlivé prac	coviště:					
	S Licence server					<u>a</u>	×
	Database						
	IP address / hos	t name: 127.	0.0.1				
	Network port:	5928	35				
	🗆 Use as licence	e server					
	Product:	Number of floati	Software version	Licence number	Running time:		
	Add			accept n	clo	se n	
						-< ¹ / ₃)-	

seca 115 Síť: Zadejte IP-adresu + port serveru v síti seca 115,

viz příručka administrátora

Licence server					a ×
atabase					
address / hos	t name: [192.	168.2.143			
etwork port:	5925	85			
Use as licence	e server				
Product:	Number of floati	Software version	Licence number	Dumning time:	
	Wantber of Hoad	Jonwale version	Licence number	nanning une.	
	rander of ridat			n onning dire.	

2. Aktivace licence jednotlivého pracoviště

Krok č. 1

Spuštění softwaru seca



Krok č. **2**

Zadání instalačního kódu (viz licenční nálepka na vnitřní straně obalu DVD) a generování registračního kódu





Zadání generovaného registračního kódu na internetové adrese seca pro obdržení aktivačního kódu

	Licence server	
	20 use node-locked licence	
	Installation code veptildbu eb07c59aus5wallar	
	Your test license expires in 23 days	
	Generate registration code	
	Registration code V060081 17X0F20TF8045X205A/4E	
	http://www.inter.com/actional activates the license. You can also upter by phone	
	Activation code	
	use floating licence	
	OK. Cancel	
■ brick brane (infolders the gain bodient, aucht () Jonak = () = (2) Adapte (2) Han — max Google = ()	anne banne bann Bender (an) Sel () han gerenne () () + 2 = + 2 3 mar () Sel () han gerenne () () () () () () () () () () () () ()	
antita interaction and a second seco	Maria de la constante de la c	AND B
Santhanan Aranan Santhanan Jarin () Jana - () - () Mang (r. 1990	Contraction	
Excellation of 20 collect per la policie (per la policie) design [c traps.com and Googles [Contraction C	
Para hann an dan sa Da bahan pata Ange (2 mpine pata Congr∈	And A land to find the set of the	
Den and an		
Den and an		
Balance of an original and an		

Krok č. 4

Obdržení aktivačního kódu prostřednictvím e-mailu, zadání a ukončení aktivace licence

Licence server	
	I use node-locked licence
Installation code	vgdtoxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx
	Your test license expires in 23 days.
	Generate registration code
Registration code	VCD08XXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX
	Please enter the license key online at <u>http://register.seco.com</u> in order to activate the license. You can also register by phone.
Activation code	x658UX000000000000000000000000000000000000
	Use floating licence

3. Změna jazyka dialogu



Spuštění softwaru seca





Vyvolejte dialogové okno "Settings"





Konfigurace bezdrátového provozu Stručný návod

Konfigurace bezdrátové sítě je nejrychlejším řešením pro seřízení přenosu dat mezi zařízeními seca 360° wireless. **Rada:** Pro rychlý přenos dat doporučujeme připojení mBCA k počítači prostřednictvím Ethernetu. Informace viz příručky administrátora softwaru a mBCA.

Krok č. 1

Připojení USB-bezdrátového adaptéru seca 456 a spuštění počítače



Krok č. **2**

Spuštění softwaru seca

Maintenance Microsoft SQL Server 2008	Music
Microsoft SQL Server 2008 R2 Microsoft Windows Performance Toolki	E Computer
Microsoft Windows SDK v7.0 Microsoft Windows SDK v7.1	Control Panel
seca software Cls configuration	Devices and Printers
S seca analytics 115	Default Programs
VMware V	- Help and Support
🕑 (2) 🔁 💽	



Vyvolejte dialogové okno "Wireless configuration" (Konfigurace bezdrátové komunikace)





Vyhledávání bezdrátových kanálů (automaticky) a spuštění přihlášení zařízení



Krok č. 5

Uložení přihlášených zařízení a ukončení konfigurace bezdrátového provozu

zobrazi z	jte, až se apnutá zařízení.	 Vyčkejte, až se zobrazí přihlášená zařízení se sériovými čís
Wireless con	figuration	Wireless configuration
Wireless group	0	Wireless group 0
	Set up wireless group	Sat up wireless group
The following o	channels are available	The following channels are available
	Channel 1	Channel 1
	Channel 3 60	Channel 2 30 Channel 3 60
	Complete login	Complete login
Please switch o group.	on all devices for the sele 2 hireless	
		The following devices are logged in:
mBCA - Length measuring	ng device	mBCA 01515201109207 Length measuing device 01274201109207
		L
	Wireless configuration	
	Wireless configuration Wireless group 0	nifers group 1 Vinifers group 2
	Wireless configuration . Windess group 0 Detap workers group The following channels are wailable	ntiess group 1 Vinstess group 2 Det up winstess group
	Wireless configuration Wireless group 0 Betrain whiters group The following channels are available Obtained 1	ntess group 1 Verstess group 2 Bet up werters group Bet up werters group Oueneel 2 Oueneel 2
	Wireless configuration Windess group 0 Driap surversi group Channel 1 Channel 2 Ohannel 3 Ohannel 4	ministri group 1 Versiless group 2 Det us vientiest group 2 Det us vientiest group 1 Outwent 1 Outwent 2 Outwent 3 Outwent 3
	Wireless configuration Wireless group 0 Earlop wireless group Dannel 1 Channel 1 Channel 2 Congritude Tages	niese goog 1 Vieless poog 2 Det up vieriess goog 1 Det up vieries poog Onaniel 1 Deaniel 1 Deaniel 1 Onaniel 2 Deaniel 3 Deaniel 3 Deaniel 3
	Wireless configuration Wireless group 0 Bring surviews group The following charinets are mailable Outward 1 Outward 2 Outward 3 Outward 4 Outward 1 Outward 2 Description logis	niess group 1 Viesless group 2 Bit up wiesless group 2 Bit up wiesless group Orannel 2 Orannel 2 Orannel 2 Diannel 2 Diannel 2 Diannel 2 Diannel 2
	Wireless configuration Wireless group 0 Bartap minimum group The following channels are involution Obarriel 1	nfers group 1 Verefers group 2 Det up winners group 2 Det up winners group 1 Ocannel 1 Ocannel 2 Biter top markets group 2 Ocannel 2 Det up winners group 1 Ocannel 2 Det up winners group 1
	Wireless configuration Wireless group 0 Brisp workers group Outmost 1 Complete large Complete large The following devices are legged in: mBCA Lange devices 01555051000007	ministi group 1 Vinisis group 2 Bit quiveriest group 2 Dianel 2 Dianel 2 Dianel 2
	Wireless configuration Wireless group 0 Brisp wareless group The following chartels are inslabilit Complete login Complete login The following devices are ingged in: Brisp Amount of Station control March are ingged in: Brisp Amount of Station control	nina ipoo 1 Vinisa pop 2 Datas kinisa pog Datas kinisa pog Datas ipoo 1 Datas 1 Datas 1 Datas 2 Datas 1 Datas 2 Datas 1 Datas 2 Datas 1
	Wireless configuration Writess group 0 Befax invites group Obtained 1 Obtained 2 Obtained 3 Obtained 4 Obtained 5 Obtained 5 Obtained 6 Obtained 6 Obtained 7 Obtained 7 Obtained 8 Obtained 9 Obtained 9 <td>verse group 1 Verses group 2 Den erse group 1 Uerses group 2 Den erse group 1 Den erses project Den erse 1 Den erses project Den erses project</td>	verse group 1 Verses group 2 Den erse group 1 Uerses group 2 Den erse group 1 Den erses project Den erse 1 Den erses project Den erses project
	Wireless configuration Wireless group 0 Bet up workerst group The following channels are walked (channels 1) Complete logie The following devices are logied in: March are searing device Stateward 1 Complete logie The following devices are logied in: March are searing device	minist group 1 Vindess group 2 Tert as wardess group 2 Tert as wardess group 1 Outward 1 Outward 2 Outward 2 Outward 2 Outward 2 Description
	Wireless configuration Wireless group 0 Betap wareless group The belowing chartels are inslatele Complete logar Complete logar The belowing devices are injeged in: BEAD Complete logar The belowing devices are injeged in: BEAD Dis2601100007 BISED Dis2601100007	

Vyvolání uživatelské dokumentace

	S seca ar	alytics 115	j			
	File Ed	lit Extras	R.			
seca analytics 11	5		S seca	analytics	115	
File Edit Extras	? Product inform Instructions fo Administration	nation r use n manual	File	Edit Extr	as ?	Product information Instructions for use Administration manua



- Zvolení tiskárny pro mBCA
-

seca analytics 115

seca analytics 115

Návod k použití

- Správa spisů pacientů seca
- Provádění měření
- Vyhodnocení výsledků vyšetření
- Tisk zpráv vyšetření
- Odeslání seca spisů pacienta k seca mBCA
-

Medical Measuring Systems and Scales since 1840

seca gmbh & co. kg Hammer Steindamm 9–25 22089 Hamburg • Germany phone +49 40 20 00 00 0 fax +49 40 20 00 00 50 technicalservice.de@seca.com

All contact data under www.seca.com



Mexico, Austria, Poland, Middle East, Brazil and Finland

Seca[®] Precision for health