seca analytics **115** Stručný návod na inštaláciu



Stručný návod (seca analytics 115





Stručný návod na inštaláciu Zriadenie pracoviska

Tento stručný návod vám ukáže, ako možno rýchlo a jednoducho uviesť softvér seca analytics 115 spolu s meracími prístrojmi seca 360° wireless do prevádzky.

Potrebné sú na to len 4 kroky:

1. Inštalácia softvéru	s. 3 – 5
2. Schválenie licencie pre jeden počítač	s. 6 – 8
3. Zmena jazyka dialógových okien	s. 9 – 10
4. Konfigurácia rádiového spojenia	s. 11 – 14

Pokyny na vykonanie meraní, vyhodnocovanie výsledkov vyšetrení, prevádzku siete a administráciu používateľov, ako aj na všetky ďalšie funkcie softvéru, nájdete v podrobnej dokumentácii pre používateľov, ktorá pozostáva z príručky pre administrátora a návodu na používanie na DVD. Aj k týmto dokumentom vám tento stručný návod ukáže cestu. **Na inštalá**ciu potrebujete na počítači práva administrátora.

1. Inštalácia softvéru

Krok č. 1

Vložte DVD. **Informácia:** Ak sa DVD automaticky nespustí, otvorte adresár DVD v programe Windows Explorer a dvojitým kliknutím kliknite na "Setup" (Inštalácia).

😰 seca analytics 115 Setup	
Seca	seca analytics 115 Welcome Setup will install seca analytics 115 on your computer. Click install to continue, options to set the install directory or Close to exit. Select option: Paralytics 115: complete Paralytics 115: complete Paralytics 115: client Seca analytics 115: liente seca analytics 115: liente seca analytics 115: liente seca analytics 115: liente Paralytics 115: liente seca analytics 115: lient
	Ļ
🖟 seca analytics 115 Setup	
SECA	seca analytics 115 Welcome Setup will install seca analytics 115 on your computer. Click install to continue, options to set the install directory or Close to exit. Select option: @ analytics 115: complete @ analytics 115: complete @ analytics 115: complete @ analytics 115: server seca analytics 115 license terms. Version 1.4.0.0 @ garee to the license terms and conditions: @ options @ pastall @ close



Sledujte pokyny na obrazovke.





Po novom spuste	ení systému sa auto	omaticky zobrazí:				
Je	edno pracovisko:					
S	Licence server				a ×	
	Database Riaddross / bost namo:	kazaar				
N	Vetwork port:	59285				
	Use as licence server					
	Product: Number	of floati Software version	Licence number	Running time:	_	
(Add		accept	close		J

sieť seca 115: IP adresu + port servera zadajte do siete seca 115, pozrite príručku administrátora

S Licence server					e ×
Database IP address / hos Network port: □ Use as licence	t name: 192. 5928 e server	168.2.143			
Product:	Number of floati	Software version	Licence number	Running time:	

2. Schválenie licencie pre jeden počítač

Krok č. 1

Spustite softvér seca.



Krok č. **2**

Zadajte inštalačný kód (pozri nálepku s licenciou na vnútornej strane obalu DVD) a vygenerujte registračný kód.





Na internetovej stránke spoločnosti seca zadajte vygenerovaný registračný kód, aby ste získali aktivačný kód.

	Licence server		
	St use and a bolto finance		
	Installation orde useful addr. sht07-50au-5us0-		
	Your test license opies in 23 days		
	Generate registration code		
	Registration code V060001 7x0F20TF80/H5x205A74E		
	Please enter the Scense key online at http://www.incense.com.in.order to activate the		
	icense. You can also here by phone		
	Activation code		
	use floating licence		
	OK.	Cancel	
	+		
2 tents bronne dickshown gane tracketer gante () brais - () = () straps - () = () straps - () = ()	Constant Sectors Control Cont : E A 2 Such and Provide (2) (2 + 2) (2 + 1) (2) (3) Sectors 2 + (2) (2) (2) (2) (2) (2) (2) (2) (2) (2)	n - Canadan Inte ⁻ O See	2015. 20 19 E C
€ Bata barra d'Arten Gra Batana gent Quanta (c) / n a Maga (c) ang (c) Grages (The second secon		
€ matajana (no kon ma Gra bana angela Angela (° man-sen Grages (And Anton (power (power (power (power (power))) (power (power)) (power)	s - Ourstein v (ore * @ See s - Courses (Seed tropage _ 3	
2 Text & sense that determine their behavior. Just C Antego (2) Text - C → C → C C → (2) - C → C → C → C → C C → (2) - C → C → C → C → C → C → C → C → C → C	And A constrained a constraine	a - c∑kanata (une * @ Ver m - c∑kanata Senetinguage _ 2 D	2010
2 Test & some data datase Gene (bedress), specific Allenge (2), specific (2), specific (2), specific Congetes ▼	And the second s	a - c∑kanata (une * 0 tra se - c∑kanata (Seectorgaage _ 2 D	2
Provi to some data data yang Grave Bordenia, yang Alanga (2), stratu - Carlos (2), stratu Grav gita > []		B Stratucki (or * * * * * * * * * * * * * * * * * * *	2
Provi haven dar denne Gene (bedreiser, genet) (Ausgelig) (⊂ transmission Genetics)		a - Cyberger Seectarysey ≥ D	2
■ Text & some data data are Units: Boards - Saint - Saint Artiging (2): Para - Saint - Saint Congets =		I - Carrier Marken (1997 - 19 Ver I - Carrier Marken (1997 - 19 Ver Seectargaage - 2 Di	2
■ Text & sense that determs (See (b) between (seek)) → A single (c) = type Grouple × [n - Cherten (or * 0 to Sectoryage 2	2
2 Text 3 men dar datars (San (bolden), jackt) Alappi (2, men Grogice)		n - Cherten (ov * 0 ve	2

Krok č. 4

Vložte aktivačný kód, ktorý ste dostali e-mailom, a ukončite schválenie licencie.

Licence server		
	🗷 use node-locked licence	
Installation code	vgdtoxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxxx	
	Your test license expires in 23 days.	
	Generate registration code	
Registration code	VCDOBXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXXX	
	Please enter the license key online at <u>http://registerseca.com</u> in order to activate the license. You can also register by phone.	L
Activation code	x658UX000000000000000000000000000000000000	
	U use floating licence	
		-

3. Zmena jazyka dialógových okien

Krok č. 1

Spustite softvér seca.





Otvorte dialógové okno "Settings" (Nastavenia).



Konfigurácia rádiového spojenia Stručný návod

Konfigurácia rádiovej siete je najrýchlejšia cesta, ako vytvoriť prenos údajov medzi prístrojmi seca 360° wireless. **Tip:** Ak chcete ešte rýchlejší prenos údajov, odporúčame vám pripojiť zariadenie mBCA k počítaču cez ethernet. Potrebné informácie nájdete v príručkách pre administrátora softvéru a prístroja mBCA.

Krok č. 1

Pripojte USB rádiový adaptér seca 456 a spustite počítač.



Krok č. **2**

Spustite softvér seca.

IronPython 2./		
Microsoft SQL Server 2008	Music	
Microsoft SQL Server 2008 R2	Computer	
🥼 Microsoft Windows Performance Toolki 🗉		
Microsoft Windows SDK v7.0	Control Panel	
seca software	Devices and Printers	
Cls configuration	benees and minters	
S seca analytics 115	Default Programs	
Startup	Help and Support	
🔰 VMware 🖬 👻		
1 Back		
Search programs and files	Shut down	



Otvorte dialógové okno "Wireless configuration" (Konfigurácia bezdrôtového spojenia).





Vyhľadajte rádiové kanály (automaticky) a spustite prihlásenie prístroja.





Prihlásené prístroje uložte a ukončite konfiguráciu rádiového spojenia.

prihláser	ijte, kým sa zobrazia é prístroje.	 Počkajte, kým sa zobrazia prihlásené prístroje so sériovými číslami.
Wireless cor	ofiguration	Wireless configuration
Wireless group	p 0	Wireless group 0
	Set up wireless group	Bet up wireless group
The following	channels are available	The following channels are available
	Channel 1	Channel 1
	Channel 3 60	Channel 2 30
	Complete login	Complete lógin
Please switch	on all devices for the select	
Shoah.		The following devices are logged in:
mBCA Length measuri	na device	mBCA 01515201109207
		Configuration of the contract
	bits down and have down	
	Wireless configuration Wireless group 0 Wireless gr	top 1 Winites group 2:
	Wireless configuration Windexs group 0 Worless pro Defuguration	No 1 Vinites prop 2.
	Wireless configuration Wireless group 0 Wireless group Barlog workers group	top 1 Virelites prop 2 Bet up weekers prop
	Wireless configuration Wireless prop 0 Wireless pro Berlig-binderse prop The following channels are available Channel 1 Channel 2 Comment 2 Channel 2 Channe	000 1 Vireless proup 2. Bertug werkess proup Ouernel 1 Dannel 1 Ouernel 2 Ouernel 2
	Wireless configuration Windexs group 0 Wireless group Barlug minimizers group	top 1 Winkes group 2 Bet deveniers group Object 1 Object 1 Object 2 Object 1 Object 2 Object 1 Object 2 Object 2 Object 2 Objec
	Wireless configuration Wireless group 0 Werless group Berlap waterer group	top 1 Virsières proup 2 Brit qui versières proup Osameil 1 Danneil 1 Osameil 2 Osameil 2 Osameil 2 Danneil 2 Danneil 2 Danneil 2 Danneil 2
	Wireless configuration Wireless group 0 Wireless pro Britigs wireless group Image: Constraint of the Constraint	top 1 Virsiers prop 2 Bet geveniers prop Ouronel 1 Ouronel 1 Ouronel 2 Dannel 2 Dannel 2 Bartigen Burtigen
	Wireless configuration Windexs group 0 Wireless group Brisp winnings group Image: State	top 1 Versies group 2 Det up werkers group 2 Octarent 1 Octarent 1 Octarent 1 Octarent 2 Octarent 2 Octarent 3 Baar tops
	Wireless configuration Wireless group 0 Werkess group Berlap warders group Image: Comparison of the Inflowing device of the Inflowing device are higherd in: Complete high Image: Complete higherd in: Darword 3 Image: Complete higherd in:	top 1 Virsières prove 2 Bet qui versières prove Outroit 2 Datamei 1 Outroit 2 Outroit 2 Outroit 2 Outroit 2 Outroit 2 Outroit 2 Outroit 2 Outroit 2 Outroit 2 Outroit 2 Outroit 2
	Wireless centiguration Wireless group 0 Wereless group Betage wireless group Image: Comparison of the State of the	top 1 Virsies group 2 Det us winners group Ouenes 1 Ouenes 1 Ouenes 2 Baarlingen Der lagen
	Wireless cestiguration Wrotess group 0 Wrotess group Britig builders group Image: 1 Dearent 2 Image: 1 Dearent 3 Image: 1 Consulter traps Image: 1 Dearent 3 Image: 1 Dearent 4 Image: 1 Dearent 5 Image: 1 Dearent 7	top 1 Versies group 2 Det up versies group Outroni 1 Det up versies group Outroni 2 Outroni 1 Outroni 2
	Wireless configuration Wireless group Berlag hearders or section Dataset Dataset Complete large Complete large	top 1 Virsites prop 2 Birl quiverieur prop Ocarent 1 Ocarent 1 Ocarent 2 Ocarent 2
	Wireless cestiguration Weekes group Bridge winniess group Weekes group The following channels are weishable Dearnel 1 Conserved 1 Dearnel 3 Conserved 1 Dearnel 4 Dearnel 3 Dearnel 4 Conserved 1 Dearnel 4 Dearnel 3 Dearnel 4 Dearnel 4 Dearnel 4 Dearnel 5 Dearnel 4 Dearnel 5 Dearnel 4 Dearnel 7 Dearnel 4 Dearnel 7 Dearnel 4 Dearnel 7 Dearnel 4 Dearnel 7 Dearnel 4 Dearnel 8 Dearnel 4 Dearnel 9	top 1 Versies poop 2 Tet querenters prope Canceral 2 C

Spustenie dokumentácie pre používateľov

	S seca analy	tics 115	ī			
	File Edit	Extras	R.			
seca analytics 115 File Edit Extras	?		S seca File	analytics Edit Ext	115 ras ?	
	Product information Instructions for use Administration man	ual				Product information Instructions for use Administration manua



- Voľba tlačiarne pre prístroj mBCA
-

seca analytics 115

seca analytics 115

Návod na používanie

- Administrácia záznamov pacientov seca
- Uskutočňovanie meraní
- Vyhodnocovanie výsledkov vyšetrení
- Tlač správ s výsledkami vyšetrení
- Odoslanie záznamu pacienta seca do prístroja seca mBCA
-

Medical Measuring Systems and Scales since 1840

seca gmbh & co. kg Hammer Steindamm 9–25 22089 Hamburg • Nemecko Telefón +49 40 20 00 00 0 Fax +49 40 20 00 00 50 technicalservice.de@seca.com

All contact data under www.seca.com



Mexico, Austria, Poland, Middle East, Brazil and Finland

Seca Precision for health