

seca EMR Module for the Windows Clipboard

The Clipboard seca EMR module enables you to copy measurements to Windows' Clipboard. Each time the *send to EMR* button is pressed, a formatted string containing the weight and/or height is copied to the Windows Clipboard. Then you may paste this string into an application of your choice.

Configuration

The configuration screen for the module looks like this:



You may specify different formats for these cases:

- There is only a weight

- There is only a height
- Both weight and height are present

In each of these formats you may use placeholders for measurements and date/time:

- {0} denotes the weight
- {1} denotes the height
- {2} denotes the date/time

Finally, using the *Append newline* setting, you can specify whether or not a newline should be appended to the end of the string.