

Medical Measuring Systems and Scales for

# HOSPITAL CARE



seca weighing and measuring solutions are designed to endure the demanding conditions of daily medical work life with ease and provide utmost accuracy even after years of service life. With seca, you can be sure that medical staff will always experience reliable state-of-the-art technology that features highest functionality, patient safety and future-proof features such as EMR-validation.



seca 676

EMR-validated wheelchair scale with handrail

- + Large wheelchair accessible platform
- + Convenient and foldable to save space
- + Mobile due to transport wheels
- + Handrail as a standing aid
- + Non-slip rubber mat
- + Two ramps included




Capacity	360 kg / 800 lbs
Graduation	50 g / 0.1 lb
Functions	TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, Auto-CLEAR, kg/lbs switch-over, damping
Accessories	+ Interface module seca 452



seca 954

EMR-validated chair scale with precise graduation

- + 50 g graduation for exact measurement while sitting
- + Fold-up arm and footrests for easy access
- + Sophisticated engineering and robust construction for a long lifespan
- + Suitable for patients affected by obesity with up to 300 kg weighing capacity



Capacity	300 kg / 660 lbs
Graduation	50 g < 150 kg > 100 g 0.1 lb < 330 lbs > 0.2 lb
Functions	TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, kg/lbs switch-over, damping, automatic switch-off
Accessories	+ Interface module seca 452 + Battery for longer mobile use seca 454

## seca 286



### EMR-validated ultrasonic measuring station with voice guidance

- + Voice-guided self-measurement assists personnel
- + Outstanding precision and high-speed measuring thanks to perfected ultrasound technology
- + Precise 50 g graduation at a very high capacity of up to 300 kg
- + High-quality workmanship and modern design with safety glass and an aluminum profile



#### Scale

Capacity	300 kg / 660 lbs
Graduation	50 g / 0.1 lb
Functions	Pre-TARE, mother+child function, HOLD, Auto-HOLD, Auto-BMI, Auto-CLEAR, CLEAR, kg/lbs switch-over, cm/inch switch-over, damping, SEND
Accessories	+ Interface module seca 452 + Animal sticker seca 487

#### Ultrasonic height measurement

Measuring range	60–210 cm / 1 ft 12"–6 ft 10"
Graduation	1 mm / 1/8"
Function	Automatic calibration

## seca 787

### EMR-validated column scale with eye-level display

- + Medical-grade, robust scale for years of use with high patient volume
- + Digital measuring rod with an extra long measuring range
- + Two-line display shows height and weight
- + Optional: Handrail seca 477 provides extra patient security
- + Optional: Panda bear facilitates measuring and weighing of children

Image contains barcode scanner from a third-party vendor.



#### Scale

Capacity	250 kg / 550 lbs
Graduation	100 g / 0.2 lb
Functions	TARE, HOLD, Auto-HOLD, damping, kg/lbs switch-over lock, cm/inch switch-over, Auto-BMI, automatic switch-off
Accessories	+ Interface module seca 452 + Handrail seca 477 + Panda seca 459 + Power adapter 401

#### Measuring rod (included)

Measuring range (total)	10–230 cm / 4–90"
Measuring range (digital)	122–230 cm / 48–90"
Graduation	1 mm / 1/8"

## seca 777

### Digital column scale with eye-level display

- + Medical-grade, robust scale designed for years of use with high patient volume
- + Built-in measuring rod with 10 – 230 cm measuring range
- + USB interface for EMR integration solutions or PC connectivity
- + Large rubber coated wheels for excellent movability
- + Optional: Handrail seca 477 provides extra patient security
- + Optional: Panda bear facilitates measuring and weighing of children

Image contains optional handrail seca 477.



#### Scale

Capacity	250 kg / 550 lbs
Graduation	100 g / 0.2 lb
Functions	TARE, HOLD, kg/lbs switch-over lock, BMI
Accessories	+ Handrail seca 477 + Panda seca 459 + Power adapter 401

#### Measuring rod (included)

Measuring range	10–230 cm / 4–90"
Graduation	1 mm / 1/8"

## seca 876

### Digital flat scale with foot pedal

- + Precise weighing technology for medical use
- + May be used standing up thanks to the foot pedals and tap-on function
- + Support feet for increased stability
- + Can weigh 10,000 times with a single set of batteries
- + Mother+child function to weigh toddlers



Capacity	250 kg / 550 lbs
Graduation	100 g < 150 kg > 200 g 0.2 lb < 330 lbs > 0.5 lb
Functions	Mother+child function, tap-on automatic switch-on, kg/lbs switch-over, automatic switch-off
Accessories	+ Stadiometer seca 217 in combination with adapter element seca 437 + Carrying cases seca 409 / seca 415 or seca 414 (for combination with stadiometer and adapter element)