



seca mBCA 525



Medical Body Composition Analyzer for determining body composition in a lying position

- + Specially designed and validated for mobile medical use
- + Patented measuring mat: BIA in supine position is finally comfortable, precise and easy
- + Lightweight and robust touchscreen monitor allows for direct evaluation
- + Unique accuracy through validation from a whole body MRI (skeletal muscle mass) and a 4C model (fat mass)
- + Validation studies are bundled and published in the European Journal of Clinical Nutrition 2013/2017, among others

mbca.seca.com

Measurement method	8-point Bioelectrical Impedance Analysis
Measurement current	100 µA
Frequencies	9
Measured data storage	Up to 100,000 measurements
Measurement duration	30 seconds
Software	seca analytics 115
Accessories	+ Rolling stand seca 475 + Basket seca 490 0019 for stand + Carrying case seca 432
Consumables	+ Single-use electrodes 490 0020

seca 954

EMR-validated chair scale with precise graduation

- + 50 g graduation for exact measurement while sitting
- + Fold-up arm and footrests for easy access
- + Sophisticated engineering and robust construction for a long lifespan
- + Suitable for patients affected by obesity with up to 300 kg weighing capacity



Capacity	300 kg / 660 lbs
Graduation	50 g < 150 kg > 100 g 0.1 lb < 330 lbs > 0.2 lb
Functions	TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, kg/lbs switch-over, damping, automatic switch-off
Accessories	+ Interface module seca 452 + Battery for longer mobile use seca 454

Medical Measuring Systems and Scales for

RENAL CARE



seca scales are indispensable tools for monitoring fluid balance, a critical aspect of treatment. By accurately measuring pre- and post-dialysis weight, seca scales enable healthcare providers to assess fluid removal efficacy, ensuring patients' safety and well-being throughout the process. seca offers a vast range of multifunctional scales to weigh immobile patients or patients with limited mobility as precise as possible.

Image contains optional handrail seca 477.



seca 777

Digital column scale with eye-level display

- + Medical-grade, robust scale designed for years of use with high patient volume
- + Built-in measuring rod with 10 – 230 cm measuring range
- + USB interface for EMR integration solutions or PC connectivity
- + Large rubber coated wheels for excellent movability
- + Optional: Handrail seca 477 provides extra patient security
- + Optional: Panda bear facilitates measuring and weighing of children

Scale

Capacity	250 kg / 550 lbs
Graduation	100 g / 0.2 lb
Functions	TARE, HOLD, kg/lbs switch-over lock, BMI
Accessories	+ Handrail seca 477 + Panda seca 459 + Power adapter 401

Measuring rod (included)

Measuring range	10 – 230 cm / 4 – 90"
Graduation	1 mm / 1/8"



seca 684

EMR-validated multifunctional scale

- + Safe weighing of standing, sitting or wheelchair patients
- + Large, easy to access platform
- + Wheel barrier provides extra safety
- + Integrated fold down seat
- + Easy to read display
- + Mobile due to integrated transport wheels



Capacity	360 kg / 800 lbs
Graduation	50 g / 0.1 lb
Functions	TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, Auto-CLEAR, kg/lbs switch-over, damping
Accessories	+ Interface module seca 452



seca 676

EMR-validated wheelchair scale with handrail

- + Large wheelchair accessible platform
- + Convenient and foldable to save space
- + Mobile due to transport wheels
- + Handrail as a standing aid
- + Non-slip rubber mat
- + Two ramps included



Capacity	360 kg / 800 lbs
Graduation	50 g / 0.1 lb
Functions	TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, Auto-CLEAR, kg/lbs switch-over, damping
Accessories	+ Interface module seca 452



seca 656

EMR-validated electronic platform scale with innovative memory function

- + Especially spacious platform for all common wheeled systems
- + Stable construction with very high capacity
- + Transportable due to wheels and integrated handles
- + Two ramps included
- + Innovative memory function for determining precise weight of patient



Capacity	360 kg / 800 lbs
Graduation	50 g / 0.1 lb
Functions	TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, Auto-CLEAR, kg/lbs switch-over, damping
Accessories	+ Interface module seca 452 + Stand for cable remote display seca 472