

# seca 210

## Mobile measuring mat for babies and toddlers



- Easy-care robust surface, suitable for commercially available disinfectants.
- Measuring mat rolls up for space-saving storage.

## seca 210:


### Adaptable devices are preferred.

Clinics, medical practices and health care facilities have varied conditions and requirements. seca responds with a wide range of products for paediatrics – from the infantometer to measuring mat and from measuring rods to the head circumference measuring tape. seca has the right solution for every application, whether the devices are for in-patient or mobile use, in large or small spaces, every day or infrequent use. They all share the renowned seca quality which always lives up to the most demanding medical requirements.

Kind to the skin,  
washable and foldable.

The mobile measuring mat seca 210 is a light weight, spacesaving solution for easy, precise measurement of the length of babies and toddlers while lying down. The fixed head piece and the sliding foot positioner make it simple to use.

*Although rugged, the surface is easy to care for and can be cleaned with all commercially available disinfectants.*



*The measuring mat can be rolled up for space-saving storage.*

## seca 210

### Technical Data

- Measuring range: 4–39" / 10–99 cm
- Graduation: 1/4" / 5 mm
- Dimensions (WxHxD):  
49.2 x 5.5 x 11.8" /  
1,250 x 140 x 300 mm
- Dimensions convoluted (WxHxD):  
4.7 x 5.5 x 11.8" /  
120 x 140 x 300 mm
- Weight: 1.3 lbs / 575 g
- With fitting for storing on wall