# seca **899**

## 



- Separate display and operating element with infinitely adjustable tilt angle for convenient use.
- Large, easy-to-read
   3/4-inch-high LCD figures.
- BMI function to determine nutritional status.



#### seca 899:

and easy to clean.

#### What was difficult is now easy.

he seca 899 is exactly right for those looking for a scale with practical cable remote display and a handle for flexible use.

This scale has many functions that make life easier. Its simple to use



## Read the display while at your desk.

The seca 899 makes your work easy. Simply place the display on your desk and the scale in front of it. Connect the two with the 7¼ ft. cable and you're ready to work. It's also easy to read the results since the tilt angle of the display is infinitely adjustable.

## Stability even on problematic surfaces.

An especially large leveling base gives the seca 899 great stability. The display and operating element is similarly equipped with a stable rubber base and with an easy-to-use self-explanatory keyboard. Storage is also effortless, thanks to its compact dimensions and light weight.

#### Measuring with all essential functions.

The seca 899 has several convenient functions required in routine medical practice: BMI to determine nutritional status, HOLD to keep the measurement results in the display, TARE to tare the weight, e.g., to weigh a small child held in a parent's arms.





Leveling base on the bottom of the scale base for safe and secure placement.



Compact design and low weight for great carrying comfort and universal usage.



Infinitely adjustable tilt angle in display and operating element for installation on a desk or wall.



### seca **899**

#### Technical Data

- Capacity: 200 kg
- Graduation: 100 g < 150 kg > 200 g
- Dimensions (WxHxD):
   321 x 60 x 356 mm / 12.6 x 2.4 x 14"
- Dimensions/platform (WxHxD):288 x 60 x 278 mm / 11.3 x 2.4 x 10.9"
- Cable length: 2 m / 78.7"
- Weight: 4.5 kg / 9.9 lbs
- Power supply: Batteries, power adapter optional
- Functions: TARE, HOLD, BMI,
   lbs/kg switch-over, automatic switch-off
- Optional: stadiometer seca 217 combined with adapter element seca 437, switch-mode power adapter seca 400, carrying case seca 414 (for combination with stadiometer and adapter element), carrying case seca 421
- Medically approved: (III)

\_