## seca integration solutions A clear improvement to your workflow

The seca 787 and seca 797 can be integrated into your information system via the existing network. At the push of a button, measured values can be sent directly to the electronic medical record (EMR), where they are immediately available to all users.

Connected with an optional USB barcode scanner, both patient and user are identified and the measurement results are always correctly assigned. This saves time, increases process reliability and eliminates transcription errors.

Get further information on seca integration solutions at: integration.seca.com

| The seca integration workflow is supported by several system providers: |  |  |  |  |
|---|--|--|--|--|
| Epic  |  |  |  |  |
| Second Allscripts   |  |  |  |  |



| Technical data                       | Seca 777<br>TARE, HOLD, lbs/kg switch-over<br>lock, BMI | Seca 787<br>TARE, HOLD, Auto-HOLD, Damping<br>(motion filter), Ibs/kg switch-over<br>lock, automatic BMI, auto-off | SECA 797<br>TARE, HOLD, Auto-HOLD, Damping<br>(motion filter), kg/lbs switch-over<br>lock, automatic BMI, USB-port and<br>holder for scanner |  |
|--------------------------------------|---|--|--|--|
| Functions                            |   |  |  |  |
| Data transmission                    | USB interface for customer specific solution            | LAN or Wi-Fi via interface module seca 452 (optional)  | Integrated Wi-Fi module  |  |
| Power supply                         | Batteries, power adapter (optional)                     | Batteries, power adapter (optional)  | Power adapter, battery pack (optional)   |  |
| Weight                               | 12.7 kg / 28.0 lbs                                      | 12.9 kg / 28.4 lbs   | 13.1 kg / 28.9 lbs   |  |
| Measuring range (total)              | 10 – 230 cm / 4 – 90 inch                               | 10 – 230 cm / 4 – 90 inch  | 10 – 230 cm / 4 – 90 inch  |  |
| Measuring range (digital)            |   | 122 – 230 cm / 48 – 90 inch  | 122 – 230 cm / 48 – 90 inch  |  |
| Capacity                             |   | 250 kg / 550 lbs   | 250 kg / 550 lbs   |  |
| Graduation                           |   |  |  |  |
| Dimensions ( $W \times H \times D$ ) |   |  |  |  |
| Platform (W $\times$ H $\times$ D)   | 35.0 × 5.9 × 41.8 cm / 13.8 × 2.3 × 16.5 inch           |  |  |  |
| Electromagnetic radiation resistance | 60601-1-2 compliant                                     |  |  |  |

seca 777 series

## Time to say good buy!



seca headquarters +49 40 20 00 00 0 info@seca.com Hamburg, Germany seca middle east +971 564 410 292 info.ae@seca.com Dubai, United Arab Emirates seca asia pacific +60 173 711 973 info.aspac@seca.com Kuala Lumpur, Malaysia

Precision for health seca.com



## seca 777 series The price-performance champion of measuring and weighing

You don't have to compromise on precision or quality when looking for a cost-efficient scale. The seca 777 provides the perfect solution with its oversized low-profile, tilt-proof solid metal platform (no plastic) and built-in measuring rod with  $10 - 230 \text{ cm} / 4 - 90^{\circ}$  measuring range. The design sets standards for precise, quick and secure weighing and measuring. With its high capacity of 250 kg / 550 lbs, the scale is optimal for accommodating all patient types. Its large, rubber coated castors also make it convenient for staff to transport the scale when needed.

If you are looking to upgrade to a digital measuring rod or EMR integration, the seca 787 and 797 are the perfect options for you.





Features

## Low-profile solid metal platform with tilt-proof design Measuring rod (10 – 230 cm / 4 – 90 inch) Large rubber coated castors for excellent movability Display swappable to the left or right side Damping system and Auto-HOLD function Digital measuring rod (122 cm up / from 48 inch) Two-line display shows height, weight and BMI USB-port and holder for scanner EMR validated with upgradeable LAN/Wi-Fi module EMR validated with built-in Wi-Fi module





