

The first scale that knows your name





seca Scale-up Line

Developed for electronic medical records

We're not only concerned with the perfect integration of measured data into your electronic medical record, we also deliver the ideal measuring station. That way, you get a completely coordinated integration solution from one source for safer and more efficient workflows.

Absolute confidence

- + The display shows names or IDs of users and patients
- + Display-controlled data transfer
- + Voice response for patient guidance aids personnel
- + Rapid measurement process within seconds thanks to the innovative ultrasound technology
- + Elegant design featuring extreme robustness thanks to a steel handrail directly anchored to the cast iron base
- + Large platform made of seca Bearclaw Glass® virtually break-resistant glass for extreme durability and hygiene





Save time, eliminate transfer errors and achieve a zeroincident environment. Transfer height and weight directly into the Electronic Medical Records (EMR) and revolutionize your workflows. seca scales and measuring stations can be integrated into all EMR systems. Feel free to contact us.





















Always the highest quality equipment for your needs



EMR-validated handrail

scale with ID-Display





- + flat scale seca 651
- + Approval class: (III)

Order code: SECA 651

seca headquarters

22089 Hamburg

Germany

seca gmbh & co. kg

Hammer Steindamm 3-25

- + Capacity: 300 kg
- + Graduation: 50 g < 150 kg > 100 g
- + handrail seca 455

+ flat scale seca 655


- + Approval class: (II)
- + Capacity: 300 kg
- + Graduation: 50 g < 150 kg > 100 g

+ Capacity: 300 kg

- + Graduation: 50 g < 150 kg > 100 g
- + Measuring range: 60 220 cm
- Graduation of measuring range: 1 mm

Order code: ON SMM AUT NN





+ flat scale seca 655

+ Approval class: (II)

- + ultrasonic height measurement seca 257
- + flat scale seca 655 + handrail seca 455
 - ultrasonic height measurement seca 257

EMR-validated ultrasonic

measuring station with

ID-Display and handrail

- + Approval class: (III)
- + Capacity: 300 kg
- + Graduation: 50 g < 150 kg > 100 g
- + Measuring range: 100 220 cm
- + Graduation of measuring range: 1 mm

Order code: ON SHM AUT NN

Order code: ON SHM ANT NN

seca united kingdom seca ltd.

40 Barn Street Birmingham B5 5QB England

headquarters

+49 40 20 00 00 0

united kingdom

+44 121 64 39 34 9

seca.com/scaleup