

Minimum PC Hardware and Software Requirements

- Microsoft-Betriebssystem®: Windows® 7, Windows® 8, Windows® 10, Windows® 11, Windows® Server 2008/R2 / Windows® Server 2012/R2, Server® 2016, Windows® Server 2019
- Processor: 1.2 GHz or higher
- Free hard disk space: min. 1 GB
- Free RAM: min. 512 MB RAM
- Connections for seca devices: USB 2.0 or higher, RJ 45 (LAN) port or adapter
- Monitor: 1024x768, High Color (16 Bit), 32 Bit
- Program for unzipping zip archives
- Program for viewing PDF files (PDF version 1.6 or higher)
- .Net Framework 4.0 (autoinstall)

NOTE:

Even if the above minimum HW requirements have been validated, we recommend a PC with an Intel Core i5 processor 2.7GHz or higher, 8GB RAM or higher.

NOTE:

Your mBCA device will be connected to the PC with a LAN cable. You need a free LAN port on your PC. If you don't have a free LAN port, or your PC doesn't have a LAN port, you'll need to use a USB-LAN adapter. This is not included in the scope of delivery and should be procured separately.

NOTE:

The seca 115 analytics 115 software is not compatible with seca 101 EMR flash or seca 103 software installation on the same PC/server