#### seca 799

### Electronic column scale with BMI function

- + Particularly robust and compact column scale
- + Includes BMI function, TARE and HOLD
- + Robust weighing platform made of steel
- + Easily movable with transport wheels
- + Can be powered with batteries for mobile use or a power adapter for stationary use
- + Expandable into a complete measuring station with an optional height measuring rod



Image contains optional measuring rod seca 220.

Scale	
Approval class	
Capacity	200 kg
Graduation	100 g < 150 kg > 200 g
Functions	TARE, HOLD, BMI, weighing range switch-over, automatic switch-off
Accessories	+ Measuring rod seca 220 + Power adapter 447 + Switch-mode power adapter 400

#### Measuring rod seca 220 (optional)

Measuring range	60-200 cm
Graduation	1 mm

seca international seca gmbh & co. kg Hammer Steindamm 3–25 22089 Hamburg Germany

T +49 40 20 00 00 0 F +49 40 20 00 00 50 E info@seca.com



#### seca 899

Г

L

### Digital flat scale with cable remote display

- + Cable remote display for wall mounting or table installation
- + Robust and durable seca quality
- + Support feet for strong stability
- + Can be powered with batteries for mobile use or a power adapter for stationary use



Approval class	
Capacity	200 kg
Graduation	100 g < 150 kg > 200 g
Functions	TARE, HOLD, BMI, automatic switch-off
Accessories	<ul> <li>+ Stadiometer seca 217 in combination with adapter element seca 437</li> <li>+ Power adapter 447</li> <li>+ Switch-mode power adapter 400</li> <li>+ Carrying cases seca 409 / seca 421 or seca 414 (for combination with stadiometer and adapter element)</li> </ul>

### Medical Measuring Systems and Scales for



٦



seca.com



#### seca 956

### Chair scale for weighing while seated

- + The simple solution for patients with decreased mobility
- + Fold-up arm and footrests make sitting down easy
- + Ergonomic seat for increased patient comfort
- + Wheels include brakes to ensure the scale stays in place when needed
- + Intuitive user interface simplifies use
- + Offers all essential functions for maximum efficiency



### seca 677

## EMR-validated wheelchair scale with handrail and transport wheels

- + Large wheelchair accessible platform
- + Convenient and foldable to save space
- + Mobile due to transport wheels
- + Handrail as a standing aid
- + Non-slip rubber mat
- + Two ramps included



Approval class	
Capacity	200 kg
Graduation	100 g
Functions	TARE, HOLD, automatic switch-off



Approval class	
Capacity	300 kg
Graduation	100 g < 200 kg > 200 g
Functions	TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, Auto-CLEAR, auto weighing range switch-over, damping
Accessories	+ Interface module seca 452



### seca 675

# EMR-validated electronic platform scale

- + Generously sized platform with high capacity of up to 300 kg
- + Flexible cable remote display with 2.3 m cable length
- + Suitable for weighing persons in a wheelchair or on a chair
- + Transport wheels and integrated handles for mobile use
- + Includes ramp for easy rolling onto the scale



Approval class	
Capacity	300 kg
Graduation	100 g < 200 kg > 200 g
Functions	TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, auto weighing range switch-over, damping
Accessories	+ Interface module seca 452 + Stand for cable remote display seca 472



Image contains optional measuring rod seca 223.

### seca 645

# EMR-validated handrail scale with generous platform

- + Safe weighing due to sturdy handrail
- + Large non-slip platform with 300 kg capacity
- + Swiveling display
- + With a BMI calculator and damping system
- + Mobile due to transport wheels



Approval class	
Capacity	300 kg
Graduation	100 g < 200 kg > 200 g
Functions	TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, Auto-CLEAR, automatic switch-off, auto weighing range switch-over, damping
Accessories	+ Interface module seca 452 + Measuring rod seca 223