

seca 674

EMR-validated electronic platform scale

- + Generously sized platform with high capacity of up to 360 kg
- + Flexible cable remote display with 2.3 m cable length
- + Suitable for weighing persons in a wheelchair or on a chair
- + Transport wheels and integrated handles for mobile use



| Capacity | 360 kg / 800 lbs |
|-------------|--|
| Graduation | 50 g / 0.1 lb |
| Functions | TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, kg/lbs switch-over, damping |
| Accessories | + Interface module seca 452 + Stand for cable remote display seca 472 |

seca international seca gmbh & co. kg Hammer Steindamm 3–25 22089 Hamburg Germany

T +49 40 20 00 00 0 F +49 40 20 00 00 50 E info@seca.com



seca 952

Chair scale for weighing while seated

- + The simple solution for patients with decreased mobility
- + Fold-up arm and footrests make sitting down easy
- + Ergonomic seat for increased patient comfort
- + Wheels include brakes to ensure the scale stays in place when needed
- + Intuitive user interface simplifies use



| Capacity | 200 kg / 440 lbs |
|------------|---|
| Graduation | 100 g / 0.2 lb |
| Functions | TARE, Auto-HOLD, kg/lbs switch-over, automatic switch-off |

Medical Measuring Systems and Scales for



and precise weighing.

L

Г



seca.com

seca 777

Digital column scale with eye-level display

Image contains optional handrail seca 477.

- + Medical-grade, robust scale designed for years of use with high patient volume
- + Built-in measuring rod with 10 – 230 cm measuring range
- + USB interface for EMR integration solutions or PC connectivity
- + Large rubber coated wheels for excellent movability
- + Optional: Handrail seca 477 provides extra patient security
- + Optional: Panda bear facilitates measuring and weighing of children

Scale

| Capacity | 250 kg / 550 lbs |
|-------------|--|
| Graduation | 100 g / 0.2 lb |
| Functions | TARE, HOLD, kg/lbs switch-over lock, BMI |
| Accessories | + Handrail seca 477 + Panda seca 459 + Power adapter 401 |

Measuring rod (included)

| Measuring range | 10-230 cm / 4-90" |
|-----------------|-------------------|
| Graduation | 1 mm / 1/8" |



seca 869

Digital flat scale with cable remote display

- + Cable remote display for wall mounting or table installation
- + Robust and durable seca quality
- + Support feet for strong stability
- + Can be powered with batteries for mobile use or a power adapter for stationary use



| Capacity | 250 kg / 550 lbs |
|-------------|--|
| Graduation | - 100 g < 150 kg > 200 g 0.2 lb < 330 lbs > 0.5 lb |
| Functions | TARE, HOLD, BMI, kg/lbs switch-over, automatic switch-off |
| Accessories | + Stadiometer seca 217 in combination with adapter element seca 437 + Carrying cases seca 409 / seca 421 or seca 414 (for combination with stadiometer and adapter element) + Power adapter 447 + Switch-mode power adapter 400 |



seca 644

EMR-validated handrail scale with generous platform

- + Safe weighing due to sturdy handrail
- + Large non-slip platform with 360 kg capacity
- + Swiveling display
- + With a BMI calculator and damping system
- + Mobile due to transport wheels



| Capacity | 360 kg / 800 lbs |
|-------------|---|
| Graduation | 50 g / 0.1 lb |
| Functions | TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, Auto-CLEAR, automatic switch-off, kg/lbs switch-over, damping |
| Accessories | + Interface module seca 452 + Measuring rod seca 223 |



seca 676

EMR-validated wheelchair scale with handrail

- + Large wheelchair accessible platform
- + Convenient and foldable to save space
- + Mobile due to transport wheels
- + Handrail as a standing aid
- + Non-slip rubber mat
- + Two ramps included



| Capacity | 360 kg / 800 lbs |
|-------------|---|
| Graduation | 50 g / 0.1 lb |
| Functions | TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, Auto-CLEAR, kg/lbs switch-over, damping |
| Accessories | + Interface module seca 452 |
| | |