

seca 677

EMR-validated wheelchair scale with handrail

- + Large wheelchair accessible platform
- + Convenient and foldable to save space
- + Mobile due to transport wheels
- + Handrail as a standing aid
- + Non-slip rubber mat
- + Two ramps included



Approval class	
Capacity	300 kg
Graduation	100 g < 200 kg > 200 g
Functions	TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, Auto-CLEAR, auto weighing range switch-over, damping
Accessories	+ Interface module seca 452

seca international seca gmbh & co. kg Hammer Steindamm 3-25 22089 Hamburg Germany



seca 877

Digital flat scale with foot pedal

- + Precise weighing technology for medical use
- + May be used standing up thanks to the foot pedals and tap-on function
- + Support feet for increased stability
- + Can weigh 10,000 times with a single set of batteries
- + Mother+child function to weigh toddlers



Approval class	
Capacity	200 kg
Graduation	100 g < 150 kg > 200 g
Functions	Mother+child function, tap-on automatic switch-on, automatic switch-off
Accessories	 + Stadiometer seca 217 in combination with adapter element seca 437 + Carrying cases seca 409 / seca 415 or seca 414 (for combination with stadiometer and adapter element)

П

Medical Measuring Systems and Scales for



Г

L



seca.com



Measuring range	60-210 cm
Graduation	1 mm
Function	Automatic calibration

Measuring range	30-220 cm
Graduation	1 mm
Functions	HOLD, automatic switch-off, SEND

Approval class	
Capacity	300 kg
Graduation	50 g < 150 kg > 100 g
unctions	TARE, Pre-TARE, mother+child function, HOLD, Auto-HOLD, BMI, Auto-CLEAR, auto weighing range switch-over, damping, automatic switch-off
Accessories	 + Interface module seca 452 + Measuring rod seca 220 or panda seca 459 + Power adapter 447 + Switch-mode power adapter 400
leasuring rod sec	ca 220 (optional)
leasuring range	60-200 cm

Measuring range 60–200 cm Graduation 1 mm



seca 959

EMR-validated chair scale with precise graduation

- + 50 g graduation for exact measurement while sitting
- + Fold-up arm and footrests for easy access
- + Sophisticated engineering and robust construction for a long lifespan
- + Suitable for patients affected by obesity with up to 300 kg weighing capacity



300 kg
50 g < 150 kg > 100 g < 250 kg > 200 g
TARE, Pre-TARE, HOLD, Auto-HOLD, BMI, Auto-CLEAR, auto weighing range switch- over, damping, automatic switch-off
+ Interface module seca 452 + Battery for longer mobile use seca 454