#### seca 704

# **EMR-validated** column scale with 300 kg capacity

- + Extremely robust, highly precise seca weighing technology with 300 kg capacity
- + Non-tilting base made of heavy duty cast iron
- + Impact- and scratchresistant powder coating
- + Wide range of functions: BMI, pre-TARE, HOLD, Auto-HOLD, mother+child function, three-step damping
- + Expandable to measuring station with an optional height measuring rod

Scale



Image contains optional measuring rod seca 220.

Approval class	
Capacity	300 kg
Graduation	50 g < 150 kg > 100 g
Functions	TARE, Pre-TARE, mother+child function, HOLD, Auto-HOLD, BMI, Auto-CLEAR, auto weighing range switch-over, damping, automatic switch-off
Accessories	<ul> <li>+ Interface module seca 452</li> <li>+ Measuring rod seca 220</li> <li>or panda seca 459</li> <li>+ Power adapter 447</li> <li>+ Switch-mode power adapter 400</li> </ul>
Measuring rod see	ca 220 (optional)
Measuring range	60-200 cm

1 mm

#### seca international seca gmbh & co. kg Hammer Steindamm 3–25 22089 Hamburg Germany

Graduation

T +49 40 20 00 00 0 F +49 40 20 00 00 50 E info@seca.com



#### seca 878

Г

L

3)

Flat scale with foot switches and double display

- + Double display lets the doctor and patient read simultaneously
- + Can weigh 10,000 times with a single set of batteries
- + Mother+child function to weigh toddlers



Approval class	
Capacity	200 kg
Graduation	100 g < 150 kg > 200 g
Functions	Mother+child function, HOLD, automatic switch-off
Accessories	<ul> <li>+ Stadiometer seca 217 in combination with adapter element seca 437</li> <li>+ Power adapter 447</li> <li>+ Switch-mode power adapter 400</li> <li>+ Carrying cases seca 409 / seca 421 or seca 414 (for combination with stadiometer and adapter element)</li> </ul>

# Medical Measuring Systems and Scales for



stations showcase child-friendly designs, including the charming animal sticker. These features ensure that the weighing process is as efficient and comfortable as possible.

Г



seca.com



Image contains optional measuring rod seca 233.

Y

#### seca 376

# iF reddot design award

# Baby scale with extra large weighing tray

- + Weighing tray with high side walls for increased safety
- + Weighing and measuring in one step thanks to optional length measuring rod
- + Three-step damping for fast weighing
- + Large tray is suitable for weighing toddlers while sitting




Image contains optional measuring rod seca 232n.

### seca 336

# Electronic baby scale for space-saving use

- + Measuring and weighing can be done in a single step thanks to an optional length measuring rod
- + Flat construction for easy storage in drawers
- + Lightweight with a carrying handle for mobile use
- + Wide range of functions: TARE, BMIF, HOLD, Auto-HOLD and adjustable damping
- + Large backlit LCD display
- + Extendable with digital measuring rod seca 234 or analog measuring rod seca 232 n



Approval class	
Capacity	20 kg
Graduation	5 g < 7.5 kg > 10 g
Functions	TARE, BMIF, HOLD, Auto-HOLD, auto weighing range switch-over, damping, automatic switch-off
Accessories	<ul> <li>+ Mobile carts seca 403 and seca 402</li> <li>+ Measuring rod seca 233</li> <li>+ Power adapter 447</li> <li>+ Switch-mode power adapter 400</li> </ul>



Approval class	
Capacity	20 kg
Graduation	5 g < 10 kg > 10 g
Functions	TARE, BMIF, HOLD, Auto-HOLD, auto weighing range switch-over, damping, automatic switch-off, confirm button
Accessories	<ul> <li>+ Mobile carts seca 403 and seca 402</li> <li>+ Digital measuring rod seca 234</li> <li>+ Analog measuring rod seca 232 n</li> <li>+ Carrying case seca 428</li> <li>+ Power adapter 401</li> </ul>



### seca 416

# Infantometer for measuring babies and toddlers

- + Robust, extra large construction for precise measuring
- + Ensures correct lying position due to trough-shaped measuring surface
- + Measurement scale printed along top side for easy reading



Measuring range	33-100 cm or 33-100 cm / 13-39"
Graduation	1 mm or 1 mm / 1/16"

#### seca 285

# **EMR-validated** measuring station for height and weight

i]=

 $\mathcal{D}$ 

GOOD DESIGN

- + Measuring and weighing in one step with absolute precision
- + Precise head alignment by the seca Frankfurt line positioner
- + Precise 50 g graduation at a very high capacity of up to 300 kg
- + Wide range of height measurement: also suitablefor pediatrics
- + Three-line touchscreen for simultaneous weight, height and BMI display
- + High-quality glass platform with a heel positioner for exact three-point measurement
- + Display in the headpiece shows the height

#### Scale

Approval class	
Capacity	300 kg
Graduation	50 g < 150 kg > 100 g
Functions	Pre-TARE, mother+child function, HOLD, Auto-HOLD, Auto-BMI, TARE, CLEAR, Auto-CLEAR, auto weighing range switch-over, damping, SEND
Accessories	+ Interface module seca 452

#### Height measurement

Measuring range	30-220 cm
Graduation	1 mm
Functions	HOLD, automatic switch-off, SEND