seca 703 s

EMR-validated column scale with 300 kg capacity

- + Extremely robust, highly precise seca weighing technology with 300 kg capacity
- + Non-tilting base made of heavy duty cast iron
- + Impact- and scratchresistant powder coating
- + Wide range of functions: BMI, pre-TARE, HOLD, Auto-HOLD, mother+child function, three-step damping
- + Expandable to measuring station with an optional height measuring rod

 \mathcal{D}

S	С	а	le	

Capacity	300 kg / 660 lbs	
Graduation	50 g < 150 kg > 100 g 0.1 lb < 330 lbs > 0.2 lb	
Functions	TARE, Pre-TARE, mother+child function, HOLD, Auto-HOLD, BMI, Auto-CLEAR, kg/lbs switch-over, damping, automatic switch-off	
Included	Stadiometer	
Accessories	 + Interface module seca 452 + Panda seca 459 + Power adapter 447 + Switch-mode power adapter 400 	

Stadiometer (included)

Measuring range	6-230 cm / 3-90"
Graduation	1 mm / 1/8"

seca international seca gmbh & co. kg Hammer Steindamm 3–25 22089 Hamburg Germany

T +49 40 20 00 00 0 F +49 40 20 00 00 50 E info@seca.com



seca 874

Flat scale with foot switches and double display

- + Double display lets the doctor and patient read simultaneously
- + Can weigh 10,000 times with a single set of batteries
- + Mother+child function to weigh toddlers



Capacity	200 kg / 440 lbs
Graduation	50 g < 150 kg > 100 g 0.1 lb < 330 lbs > 0.2 lb
Functions	Mother+child function, Auto-HOLD, kg/lbs switch-over, automatic switch-off
Accessories	 + Stadiometer seca 217 in combination with adapter element seca 437 + Carrying cases seca 409 / seca 415 or seca 414 (for combination with stadiometer and adapter element)

Г			٦

L

Medical Measuring Systems and Scales for



stations showcase child-friendly designs, including the charming animal sticker. These features ensure that the weighing process is as efficient and comfortable as possible.



seca.com



seca 374

Y reddot design award

Baby scale with extra large weighing tray

- + Weighing tray with high side walls for increased safety
- + Weighing and measuring in one step thanks to optional length measuring rod
- + Three-step damping for fast weighing
- + Large tray is suitable for weighing toddlers while sitting



Capacity	20 kg / 44 lbs
Graduation	5 g < 10 kg > 10 g / 0.2 oz < 22 lbs > 0.5 oz
Power supply	Batteries, power adapter (optional)
Functions	TARE, BMIF, Auto-HOLD, kg/lbs switch-over, damping, automatic switch-off
Accessories	 + Mobile carts seca 403 and seca 402 + Measuring rod seca 233 + Power adapter 447 + Switch-mode power adapter 400



Image contains optional measuring rod seca 232.

seca 334

Mobile baby scale with flat weighing tray

- + Lightweight with a carrying handle for mobile use
- + Measuring and weighing can be done in a single step thanks to an optional length measuring rod
- + Flat construction for easy storage in drawers
- + Large, easy to read LCD display
- + Wide range of functions: TARE, BMIF, Auto-HOLD and adjustable damping



Capacity	20 kg / 44 lbs
Graduation	5 g < 10 kg > 10 g / 0.2 oz < 22 lbs > 0.5 oz
Power supply	Batteries, power adapter (optional)
Functions	TARE, BMIF, Auto-HOLD, kg/lbs switch-over, automatic switch-off
Accessories	 + Mobile carts seca 403 and seca 402 + Measuring rod seca 232 + Power adapter 447 + Switch-mode power adapter 400 + Carrying case seca 428



seca 416

Infantometer for measuring babies and toddlers

- + Robust, extra large construction for precise measuring
- + Ensures correct lying position due to trough-shaped measuring surface
- + Measurement scale printed along top side for easy reading



Measuring range	33-100 cm or 33-100 cm / 13-39"
Graduation	1 mm or 1 mm / 1/16"



seca 284

EMR-validated measuring station for height and weight

GOOD DESIGN

3)

- + Measuring and weighing in one step with absolute precision
- + Precise head alignment by the seca Frankfurt line positioner
- + Precise 50 g graduation at a very high capacity of up to 300 kg
- + Wide range of height measurement: also suitablefor pediatrics
- + Three-line touchscreen for simultaneous weight, height and BMI display
- + High-quality glass platform with a heel positioner for exact three-point measurement
- + Display in the headpiece shows the height

Scale

300 kg / 660 lbs
50 g / 0.1 lb
Pre-TARE, mother+child function, HOLD, Auto-HOLD, Auto-BMI, Auto-CLEAR, CLEAR, kg/lbs switch-over, cm/inch switch-over, damping, SEND
+ Interface module seca 452

Height measurement

Measuring range	30-220 cm / 11"-7 ft 2"
Graduation	1 mm / 1/8"
Functions	HOLD, user-defined zero setting, automatic switch-off, SEND