





Graduation

seca 216

adults

Mechanical measuring

+ Easy-to-read measurement scale

+ Choice of two units of measurement

rod for children and

+ Adjustable head positioner

+ Simple assembly



seca 206

Mechanical measuring tape

- + Space-saving, roll-up mechanism
- + Simple read-out of results through display window
- + Durable tape made of metal




Measuring range max. 138 cm / 55.3" selectable from 3.5-230 cm / 1.4-90"



Measuring range 0-220 cm or 0-87" 1 mm or 1/8" Graduation

> You can feel seca quality, it gives a good impression and pays off. The caliper gauge of a seca stadiometer runs smooth and with millimetre precision over the abrasion-resistant scale and a seca column scale gives a safe and high-quality impression, even after thousands of weighings. It is not without reason that seca is the world market leader in the field of medical weighing and measuring and thus the quality standard for every medical practice.

Medical Measuring Systems and Scales for









1 mm / 1/16"

T +49 40 20 00 00 0 F +49 40 20 00 00 50 E info@seca.com

seca.com



Image contains optional measuring rod seca 232n.



Electronic baby scale for space-saving use

- + Measuring and weighing can be done in a single step thanks to an optional length measuring rod
- + Flat construction for easy storage in drawers
- + Lightweight with a carrying handle for mobile use
- + Wide range of functions: TARE, BMIF, HOLD, Auto-HOLD and adjustable damping
- + Large backlit LCD display
- + Extendable with digital measuring rod seca 234 or analog measuring rod seca 232 n



Approval class	
Capacity	20 kg
Graduation	5 g < 10 kg > 10 g
Functions	TARE, BMIF, HOLD, Auto-HOLD, auto weighing range switch-over, damping, automatic switch-off, confirm button
Accessories	+ Mobile carts seca 403 and seca 402 + Digital measuring rod seca 234 + Analog measuring rod seca 232 n + Carrying case seca 428 + Power adapter 401



seca 384

2-in-1 mobile baby scale and flat scale for toddlers

- + Weigh babies lying down and toddlers standing up
- + Easily detach the weighing tray from the base to use as a flat scale
- + Compact, robust and lightweight: ideal for mobile use
- + High-quality baby scale with calibration approval at an economical price
- + Wide range of functions: TARE, HOLD and Breast-Milk-Intake-Function (BMIF)





Approval class	
Capacity	20 kg
Graduation	10 g < 10 kg > 20 g
Functions	TARE, BMIF, HOLD, automatic switch-off, weighing range switch-over
Accessories	 + Mobile carts seca 403 and seca 402 + Measuring board seca 417 or measuring mat seca 210 + Backpack seca 409 (for combination with measuring board or measuring mat) + Carrying case seca 414 (for combination with measuring board) + Backpack seca 431 + Carrying case seca 413

seca 799

Electronic column scale with BMI function

- + Particularly robust and compact column scale
- + Includes BMI function, TARE and HOLD
- + Robust weighing platform made of steel
- + Easily movable with transport wheels
- + Can be powered with batteries for mobile use or a power adapter for stationary use
- Expandable into a complete measuring station with an optional height measuring rod



Image contains optional measuring rod seca 220.

seca 878 dr

Flat scale with customizable label

- + Customizable label helps minimize scale misplacement
- + Inexpensive, reliable and durable: seca quality that pays off
- + Double display allows the doctor and patient to read simultaneously
- + Mother+child function for hassle-free weighing of toddlers

seca.com/doctorscale



Approval class	
Capacity	200 kg
Graduation	100 g < 150 kg > 200 g
Functions	Mother+child function, HOLD, automatic switch-off
Accessories	+ Stadiometer seca 217 in combination with adapter element seca 437 + Power adapter 447 + Switch-mode power adapter 400 + Carrying cases seca 409 / seca 415 or seca 414 (for combination with stadiometer and adapter element)

Scale

Approval class	
Capacity	200 kg
Graduation	100 g < 150 kg > 200 g
Functions	TARE, HOLD, BMI, weighing range switch-over, automatic switch-off
Accessories	+ Measuring rod seca 220 + Power adapter 447 + Switch-mode power adapter 400

Measuring rod seca 220 (optional)

Measuring range	60-200 cm
Graduation	1 mm