

Ergonomic circumference measuring tape

seca 201

- + Retractable measuring tape for precise body measurements
- + Printed scale on side makes results easy-to-read
- + Convenient one-handed operation fasten start of tape to retraction mechanism
- + Compact and easy to store



seca 201

Ergonomic circumference measuring tape



A high-quality medical device

Because circumference measuring tapes today also have to meet the highest medical requirements, we have developed the seca 201, with which the body circumference can be measured with millimeter precision. The high-quality mechanism ensures that the 205 cm long tape can be easily pulled out and snaps in precisely. The housing fits well in the hand, feels pleasant and is so robust that it can also withstand being accidentally dropped.

Measuring range	0 – 205 cm or 0 – 80 in
Graduation	1 mm or 1/8 in
Dimensions (WxHxD)	70 × 22 × 65 mm / 2.8 × 0.9 × 2.6 in
Net weight	0.05 kg / 0.1 lbs
Medical Device Regulation (MDR) certified	

seca international
seca gmbh & co. kg
Hammer Steindamm 3–25 T +49 40 20 00 00 0
22089 Hamburg F +49 40 20 00 00 50
Germany E info@seca.com

Precision for health
[seca.com](https://www.seca.com)