

Mobile stadiometer for measuring height seca 214



- ▲ Simple and easy to set up - no wall fastening necessary
- large floor plate ensures stability
- result clearly visible while measuring
- convenient and easy to transport

The convenient mobile method of measuring height

The stadiometer seca 214 is especially suitable for mobile use, e.g. for measuring children and teenagers in the course of medical examinations in schools, for use by carers and health workers on the move. But it is also ideal for doctor's surgeries and hospitals specialising in child and adolescent medicine. Furthermore, the seca 214 is suitable for measuring patients lying down. Once taken apart, the stadiometer is easy to transport anywhere.

Simple and quick to set up - no wall fastening required

The measuring rod seca 214 which can be dismantled into several pieces can be set up easily and quickly. It is mounted on the floor plate in no time at all. Further fixation, e.g. by means of a wall fastening, is no longer necessary.



Flexible positioning due to large floor plate

The stadiometer seca 214 can be used anywhere because the large floor plate provides the necessary stability. The stadiometer seca 214 is therefore the ideal piece of equipment for all those who require reliable measuring instruments at frequently changing locations.

Result clearly visible while measuring

Adjusting the smooth-running measuring slide is extremely easy due to the large head plate. A further sophisticated detail is that, as the scale is printed along the side of the measuring rod, it is easy to read off the result while measuring, thus guaranteeing precise results - up to a height of 207 cm.

Convenient and easy to transport

Practical: the stadiometer seca 214 is easy to take anywhere not only because it can be dismantled in no time at all but also because it has a practical, integrated handle.



seca 412 – carrying case for measuring rod seca 214

Technical data

seca 214

Measuring range: 20 - 207 cm / 8 - 81.5"

Graduation: 1 mm / 1/8"

Dimensions (WxHxD): 356 x 2,116 x 419 mm

Weight: 2 kg

Optional: carrying case seca 412

