

seca 354m

Digital baby scale,
also ideal as floor scale
for children

optimized



- Baby scale and floor scale for children in one.
- Weighing tray easy to detach from base.
- Breast-Milk-Intake-Function.
- Mobile and robust for transport.
- Operated by power adapter or batteries.

seca 354m:

Convertible and robust –
compact quality with dual purpose.

At home, in surgery and around the world wherever life begins: reliable weight control is decisive for a healthy start in life. The seca 354m is an all-around scale which has been developed on the basis of practical experience. It can be used not only as a baby scale but also as a floor scale for small children.



Baby scale and floor scale for children in one.

Babies can be weighed safely and comfortably in the tray of the seca 354m. The rounded shape of the tray with its slightly raised sides gives the baby a feeling of security. By detaching the tray from the base, the scale is converted into a fully functional digital floor scale with a capacity of 20 kg.

Very safe due to secure attachment of tray.

In no time at all, the baby scale can be converted into a floor scale: pressing a button unfastens the mechanism attaching the tray to the base. When locked, this secure attachment ensures that the baby scale is easy to lift and safe to transport.



The weighing tray can be removed quickly and easily at the touch of a button.

Easy to operate due to self-explanatory button functions.

The seca 354m is easy to operate due to its intuitive operating concept with self-explanatory buttons. The auto-HOLD function means that, after being weighed, the baby can first be taken care of before one has to read off the result. With the TARE function, the weight of nappies or soft mats can be reliably deducted. With the Breast-Milk-Intake-Function, one can determine how much milk the baby has drunk when breast fed. Because, at the touch of a button, the weight of the baby determined before feeding is deducted from the weight determined after feeding.



Mobile and robust for transport.

Its compact, robust construction predestines the seca 354m as the scale for all those who have to be mobile in the interest of the health of their youngest patients. The power-saving battery operation – up to 20,000 weighings with just one set of batteries – makes the scale independent of electricity power supply.



Versatile operation due to battery or power supply.

seca 354m

Technical Data

- Capacity: 20 kg / 44 lbs
- Graduation: 10 g < 10 kg > 20 g, 0.5 oz < 22 lbs > 1 oz
- Dimensions (WxHxD): 552 x 156 x 332 mm
21.7 x 6.1 x 13.1"
- Dimensions/tray (WxHxD): 525 x 80 x 250 mm
20.7 x 3.1 x 9.8"
- Dimensions/platform (WxHxD): 280 x 65 x 235 mm
11 x 2.6 x 9.3"
- Weight: 2.3 kg / 5 lbs
- Power supply: Power adapter, batteries
- Functions: TARE, BMIF, auto-HOLD, kg/lbs switch-over, automatic switch-off
- Optional: carrying case seca 414 and seca 413