

Baby scale with extra large weighing tray

seca 374

- + Weighing tray with high side walls for increased safety
- + Weighing and measuring in one step thanks to optional length measuring rod
- + Three-step damping for fast weighing
- + Large tray is suitable for weighing toddlers while sitting
- + Medical-grade seca scale technology with precise graduation



seca 374

Baby scale with extra large weighing tray



Professional scale for users, feel-good scale for babies

The understanding of babies for the necessary weight control is often kept within limits. They have just been sleeping peacefully in their parents arms and now are lying wide awake, naked and alone on a cold scale. In order to make this scenario as pleasant as possible for the littlest of little ones, the seca 374 baby scale has particularly high side walls. This provides a sense of security, protects from falling out and creates limitations. Because wide, open spaces can be frightening for babies.

Capacity	20 kg / 44 lbs
Graduation	5 g < 10 kg > 10 g / 0.2 oz < 22 lbs > 0.5 oz
Dimensions (WxHxD)	620 × 190 × 358 mm / 24.4 × 7.5 × 14.1 in
Tray (WxHxD)	615 × 130 × 275 mm / 24.2 × 5.1 × 10.8 in
Net weight	3.6 kg / 7.9 lbs
Power supply	Batteries, Power adapter (optional)
Functions	TARE, Auto-HOLD, RESET, Automatic switch-off, kg/lbs switch-over, kg/lbs switch-over lock, Adjustable damping, BMIF
Accessories	+ Analog measuring rod seca 233 + Switch-mode power adapter seca 400 + Mobile cart seca 402 + Mobile cart with drawer seca 403

Medical Device Regulation (MDR) certified



Weighing and measuring in one step

The seca 374 baby scales can be equipped with the high-quality seca 233 aluminum measuring rod. However, weighing and measuring in one work step is not the only way to save time. With the three-stage adjustable damping, the weight can be quickly determined, even with a lively baby, by filtering out the movements. In addition, the baby can leave the scale before the value has been recorded, because the HOLD function records the result.



The details make the difference

The seca 374 proves itself with many details that make it ideal for use in medical settings. For example, the seca 374 has rubberized feet for a non-slip setup. The ergonomic weighing tray is cast as a single piece without recesses or grooves and is easily disinfected. The display is positioned so that it can be read without distracting attention from the infant.