

Light, space-saving and
stable measuring board
also ideal for mobile use

seca 417

- + Simple operation for quick mobile use
- + Folds up for space-saving storage
- + High-quality materials for long service life
- + Smooth surfaces for easy cleaning



seca 417

Light, space-saving and stable measuring board also ideal for mobile use



Easy to read

The seca 417 stadiometer is ideal for mobile use to check the nutritional and developmental status of babies and toddlers, e.g., during serial examinations. It is precisely here that it is important that measuring results can be read quickly and easily. This is guaranteed by red reading lines and high-quality scale printing, which does not fade even after years of intensive use.

Easy to operate

The two-part seca 417 measuring board folds up via a hinge. It has an integrated head positioner and a separate foot positioner, which can be moved easily and precisely on a gently raised side rail. For measuring purposes, it is sufficient to align the child's head with the head positioner, carefully stretch the legs and guide the foot positioner to the feet.



Quick cleaning

The rounded edges and smooth surface are for baby's safety and allow for easy cleaning and disinfecting. The foot positioner is detachable so that pathogens cannot collect in the side rails. Commercially available cleaning agents and disinfectants cannot affect the material made of skin-friendly ABS plastic and the printing on the scale.

Measuring range	10 – 100 cm / 4 – 39 in
Graduation	1 mm / 1/16 in
Dimensions (WxHxD)	1,109 × 115 × 333 mm / 43.7 × 4.5 × 13.1 in
Dimensions when folded (WxHxD)	577 × 115 × 333 mm / 22.7 × 4.5 × 13.1 in
Net weight	1.6 kg / 3.5 lbs
Accessories	+ Backpack seca 409 + Carrying case seca 412 + Carrying case seca 414

Medical Device Regulation (MDR) certified

seca international
seca gmbh & co. kg
Hammer Steindamm 3–25 T +49 40 20 00 00 0
22089 Hamburg F +49 40 20 00 00 50
Germany E info@seca.com

Precision for health

[seca.com](https://www.seca.com)