
seca 472 Mobile stand for cable remote displays of seca scales.



- Stable due to non-slip rubber feet.
- Sturdy, durable construction.



seca 472:

The stable solution for cable remote displays.

his one is not so easy to bowl over: the stand seca 472 is the transportable solution for holding the display and control element of the digital scales seca 634, 674 and 656.

Firm hold and easy to transport.

The stand seca 472 is not only a practical alternative to the table support or wall fixture. Equipped with non-slip rubber feet, it also replaces the former model on transport castors and can be positioned everywhere quickly and securely.

Sturdy, durable construction.

The stable stand is made of high grade steel tubing and is thus of durable quality. The display and operating control element is held firmly by a robust clamp. The height and inclination angle of the display element can be adjusted to individual requirements. It can therefore be operated comfortably both from a seated or standing position.








seca 472

Technical data

- Dimensions (WxHxD): 18.1 x 40.7 x 18.1 inch / 460 x 1035 x 460 cm
- Weight: 17.2 lbs / 7.8 kg





Г