

seca 674

High capacity digital platform scale with transport castors



- Drive-up ramp seca 470 included.
- Suitable for weighing patients in a wheelchair or chair.
- Flexible cable remote display.
- Transport castors.

seca 674:

Weighing without restrictions.

Special situations demand special scales. Not every patient can step onto a scale quickly and easily. Patients who are overweight and those who have limited mobility will need the help of a convenient scale like the seca 674. The platform scale serves a wide range of patients and eases the routine work of doctors and nurses.

Large platform with high load-bearing capacity.

This generously-sized wheelchair/platform scale is multifunctional. It weighs patients in care, especially heavy patients, persons in wheelchairs and seated patients with circulatory trouble, even during dialysis. The flexible cable remote display can be placed anywhere (e.g. on a desk or table) for convenient read-out and operation of the controls. Wall mount and seca 470 drive-up ramp is included. Transport castors and a handle make the seca 674 mobile. The Pre-TARE function simply subtracts the added weight of the chair or wheelchair and provides the weight of the patient alone.

Versatile wired remote display.

Since the display element of the seca 674 is connected to the scale by means of 2.7-meter-long cable, it can be positioned where the user can most easily read the results. An additional feature is the RS232 adapter set seca 460 (optional) which connects this scale to a printer or PC.

Optional: non-slip base.

The seca 472 stand is the stable solution for cable remote displays of seca 674. The stable construction with non-slip rubber feet replaces the former model on transport castors and can be positioned everywhere quickly and securely. The display element can easily be attached to the stand by means of a stable bracket. The height of the display is continuously adjustable up to 40 inches.

Easy to transport.

Despite its size, the seca 674 is amazingly easy to handle. When upright, the scale can be moved easily on castors and steered by means of the integrated handle. This feature makes the seca 674 suitable for mobile use.



The seca 472 stand for cable remote displays of seca scales and measuring rods.



A ramp for easy wheelchair access is included in delivery.

seca 674

Technical data

- Capacity: 800 lbs / 360 kg
- Graduation: 0.1 lbs / 50 g
- Dimensions (WxHxD):
35.6 x 2.8 x 41.9" / 904 x 70 x 1,064 mm
- Dimensions/platform (WxHxD):
29.9 x 2.2 x 36.6" / 760 x 55 x 930 mm
- Weight: 59.5 lbs / 27 kg
- Power supply: Power adapter
- Functions: TARE, pre-TARE, HOLD, auto-HOLD, BMI, auto-CLEAR, lbs/kg switch-over, damping
- Optional: Stand for cable remote display seca 472, second ramp seca 470
- One ramp included