# EMR-validated wheelchair scale with handrail and transport wheels



# seca 676

- + Large wheelchair accessible platform
- + Convenient and foldable to save space
- + Mobile due to transport wheels
- + Handrail as a standing aid
- + Non-slip rubber mat
- + Two detachable ramps included
- + EMR-validated: transmit measured data directly to an electronic medical record system or PC



# EMR-validated wheelchair scale with handrail and transport wheels



## seca 676 - wheelchair accessible weighing made easy

The seca 676 can be used as a wheelchair scale or as a "seat scale" with a chair, where its tare weight is deducted. The handrail always provides a secure feeling. It is so firmly anchored that even a patient affected by obesity can lean against it without it giving in. The measuring result remains unaffected by this. The non-slip rubber pad and the safety rail ensure additional safety, preventing the wheelchair from rolling over to the side.





### Mobile and space-saving

If the patient does not come to the seca 676, the seca 676 just comes to the patient. Folded up and with transport castors, an individual person can move the wheelchair scale to another room since the handrail serves as a handle. In order to save space, the seca 676 can be folded up and removed when required.



# seca integration solutions: How to significantly improve workflows

Integration capability means that the seca 676 is integrated into the hospital information system via the network. In this way, measured data can be sent directly to the digital patient file (EMR) at the touch of a button, where they are immediately available to all users. With an additionally connected USB barcode scanner, the patient and user are identified and the measurement results are always correctly assigned. This saves you time, increases process reliability and eliminates transmission errors. More information on seca integration solutions can be found at: integration.seca.com