seca **769**

Digital column scale with BMI function



- Site-independent as battery operated.
- Easy to transport.
- Cost-effective and low-maintenance due to low power consumption.
- seca 220 measuring rod included.



seca 769:

The intelligent solution for flexible use.

must in the medical sector: high precision and sophisticated functionality combined with cost-effectiveness. For over one hundred years, seca has continued to develop innovative scales and measuring systems to provide optimal solutions. The seca 769 is perfectly aligned with the requirements that hospitals and medical practices must achieve, is cost effective and eases the workload for doctors and nursing staff.



Intelligent functions for versatile use.

Determining the nutritional condition of the patient is essential in hospitals and doctors' practices. That is why the seca 769 has the BMI function where the height of the patient is compared with his weight, thus providing valuable information. The auto-HOLD function too reduces the everyday workload: it ensures that the result is still shown on the display after the patient has stepped down from the scale. The practitioner can therefore first attend to the patient before making a note of the weight. Additional items which are not to be included in the weight result are simply deducted by means of the TARE function.

Site-independent and easy to transport.

The seca 769 can be used wherever it is needed. As it is battery-operated, no wall outlet connection is necessary. The seca 769 is easy to move around with its smooth running transport castors which can simply be wheeled wherever it is needed.

Especially cost-effective and low-maintenance.

The seca 769 is energy-saving, thus very costeffective. The minimum power consumption of the low maintenance column scale ensures that a large number of weighings are possible with just one set of batteries. The automatic switchoff function is another reason for the low energy consumption.



All it takes is the press of a key to find out BMI.



The included seca 220 measuring rod allows simultaneous weighing and measuring.

seca **769**

Technical data

- Capacity: 450 lbs / 200 kg
- Graduation: 0.2 lbs / 100 g
- Dimensions (WxHxD):
 11.6 x 32.7 x 16.4" /
 294 x 831 x 417 mm
- Dimensions platform (WxHxD):
 10.7 x 3 x 11" /
 272 x 75 x 280 mm
- Weight: 13.7 lbs / 6.2 kg
- Power supply: batteries, power adapter optional
- Functions: TARE, HOLD, BMI, lbs/kg switch-over, automatic switch-off
- Optional: measuring rod seca 224, power adapter seca 447, switch-mode power adapter seca 400