Digital column scale with eye-level display



seca 777

- + Medical-grade, robust scale designed for years of use with high patient volume
- + Built-in measuring rod with 4 90" measuring range
- + USB interface for PC connectivity
- + Large rubber coated castors for excellent movability
- + Robust tilt-proof full metal platform, no plastic
- + 8-year warranty (4x the industry standard)
- + Optional: Handrail seca 477 provides extra patient security
- + Optional panda bear facilitates measuring and weighing of children



seca 777

Digital column scale with eye-level display



Cost-effective without any compromises in quality

The seca 777 puts unstable and inaccurate weighing and measuring to an end. It boasts a tilt-proof full metal platform and a stable single piece column. The eye-level display and ultra sturdy base provide comfort and security. It is a quality every patient will notice, as soon they sets foot on the platform. The scale is easy to move by staff due to the large rubber coated castors and a wide wheelbase.

Scale

Capacity	550 lbs / 250 kg
Graduation	0.2 lbs / 100 g
Dimensions (WxHxD)	16.1 x 53.4 x 22.7 inch / 41.0 x 135.6 x 57.6 cm
Platform (WxHxD)	13.8 x 2.3 x 16.5 inch / 35.0 x 5.9 x 41.8 cm
Measuring range	4 – 90 inch / 10 – 230 cm
Weight	28.0 lbs / 12.7 kg
Power supply	Batteries, power adapter (optional)
Data transmission	USB interface for customer specific solution
Functions	TARE, HOLD, Ibs/kg switch-over lock, BMI
Electromagnetic radiation resistance	60601-1-2 compliant
Accessories	 + Handrail seca 477 + Panda seca 459 + Power adapter 401

Handrail seca 477 (optional)

Dimensions (WxHxD)	16.9 x 44.2 x 12.2 inch / 42.9 x 112.2 x 30.9 cm
Weight	8.8 lbs / 4.0 kg



We reserve the right to make modifications.

201177704022

seca north america | east 498 Wando Park Boulevard Mount Pleasant, SC 29464 USA

T +1 800 542 7322 F +1 888 705 7397 E info.us@seca.com



Optional handrail for increased patient stability

To accommodate patients in need of extra support during the weighing process, the seca 777 is equipped with the option to add an integrated handrail (seca 477). The easy-to-mount handrail raises the level of patient security to a new level. The measuring rod is securely attached to the column and provides a wide measuring range that is suitable for pediatrics to tall individuals. In addition, we've added extra space between the column and the platform to ensure the highest comfort.



Sturdy construction and quality you can trust

As the world market leader in medical measuring and weighing, we put all our experience and knowledge in the development and design of the seca 777. This includes extensive and stringent endurance tests that simulate more than ten years of daily usage. We made sure, that even after thousands of height and weight measurements and thousands of button presses, the seca 777 will work smooth and reliable as it did on the first day.