

seca 786

Mechanical column scales with large round dial



- With hip-high column for easy read-out.
- With large black platform.
- Ready for use without electrical connection or battery.
- Integrated transport castors.

seca 786:

Classic highlight in weighing technology.

This new circular balance has properties that are tried and tested by seca over several decades but still modern today. The seca 786 is always ready for use because it needs neither electricity nor batteries. The robust scale determines weight quickly, accurately and conveniently. It is manufactured in accordance with seca top-quality standards, as expected from the world market leader in medical weighing and measuring.

Exact results, uniquely designed.

With its hip-high column, the scale adds elegance to any location. The large round dial with a red needle precisely displays the easy-to-read weight in 500 g graduation. Every patient can weigh himself effortlessly without assistance.

So high, so practical.

The seca 786 combines reading ease with space-saving dimensions. The elevated dial leaves enough room for a large platform. Its black rubber coating is not only robust but also easy-to-clean.



The separately available measuring rod seca 224 permits simultaneous weighing and measuring.

The black rubber mat on the platform is robust and easy-to-clean.



seca 786

Technical Data

- Capacity: 150 kg
- Graduation: 500 g
- Dimensions (WxHxD): 295 x 910 x 500 mm
- Weight: 6.9 kg
- Optional: measuring rod seca 224