seca **803**

Digital flat scale for individual patient use



- Two-component rubber coating for pleasant footing.
- Four stable rubber feet for a safe standing.
- 1.1 inch LCD figures for easy reading.



seca 803:

The rediscovery of substance.

irst take a look at the clear design and then at the finishing: even our simplest scale is complex in its design and robust in its construction. The seca 803 allows patients who weigh up to 330 lbs safe footing.



The rediscovery of feeling good.

Go ahead, pick up our starter scale and take a close look. You'll be impressed by its classic design and its surface gleaming with high-quality materials. The high-quality rubber coating consists of a special combination of two components which lends softness and solidity. The material can stand up to all commercially available cleaning and disinfecting products.

Professional can be so private.

Whether black or white, for professional or private use, this seca starter scale satisfies the highest demands for private use and anyone who appreciates high-quality workmanship. The seca 803 has its strengths on the outside and on the inside. Underneath the weighing platform you'll see longitudinal and cross braces set one centimeter apart, giving the scale unparalleled stability.


Lightening the load.

Thanks to its low tare weight of only 3 lbs and the handle built on the back, the scale is easy to carry and makes it suitable for users on the move. Viewing the measurements is just as easy from the unconventionally large display with 1.1 inch-high figure. No corners are cut with the seca 803.





Close ribbing for a secure foothold.





seca **803**

Technical data

- Capacity: 330 lbs / 150 kg
- Division: 0.2 lbs / 100 g
- Dimensions scale (WxHxD):
 12.4 x 1.4 x 12.8 inch
- Weight: 3.3 lbs
- Functions: Auto-HOLD, automatic switch-off
- Parts to be delivered:
 4 AAA batteries,
 4 extension feet for carpets