

Flat scale with foot switch and optional Bluetooth® interface

seca 874 +

- + Developed for mobile use and protected against environmental influences with an IP53 protection code
- + Platform with an anti-slip surface and antimicrobial plastic for optimal hygiene
- + Extremely stable steel base with a heavy-duty weight capacity
- + Convenient and easy to operate with a single foot switch for efficient use
- + LC-display can be oriented towards the patient or doctor, offering convenience and flexibility
- + Mother+child function simplifies the weighing of children in the arms of the mother or father
- + Optional Bluetooth® interface: enables error-free data transmission



seca 874 +

Flat scale with foot switch and optional Bluetooth® interface



The benchmark for mobile medical weighing

Introducing the seca 874 +, an enhanced version of our classic flat scale, setting a new standard for mobile medical use. We've enlarged the platform, optimized the foot switch for a streamlined operation, and reduced the net weight. The LC-display can be oriented based on your preferences — towards the patient or the medical staff. For enhanced functionality, an optional Bluetooth® interface facilitates error-free automatic weight reading transmission to a dedicated app.



Hygiene standards at the forefront

Incorporating an antimicrobial additive into the plastic tread surface inhibits germ spread, ensuring impeccable hygiene. The easy-to-clean surface eliminates unnecessary “dirt traps,” with rounded edges for effortless cleaning. The specially textured surface not only guarantees hygiene but also provides high slip resistance.



Robust and durable construction

Equipped with leveling feet, the seca 874 + remains stable and delivers precise results even on uneven surfaces. The steel base ensures stability while remaining lightweight for mobile use. Rigorous testing ensures durability, dust resistance, and splash water protection according IP53 protection code. The seca 874 + is ready for any challenge, even in non-air-conditioned rooms, ensuring a long-scale life.

Capacity	200 kg
Graduation	100 g
Dimensions (WxHxD)	14 x 2.4 x 15.7 in / 36 x 6.1 x 40 cm
Platform (WxHxD)	13.8 x 2.3 x 12.4 in / 35 x 5.74 x 31.5 cm
Net weight	8 lbs / 3.7 kg
Power supply	Batteries, Power adapter (optional)
Functions	TARE, HOLD, Auto-HOLD, CLEAR, RESET, Automatic switch-off, Mother+child function
Interfaces	Bluetooth® (optional)
Accessories	+ Switch-mode power adapter seca 400 + Power adapter 401

seca 874 + variants	8741009719	874320971
Capacity: 200 kg	•	•
Interface: Bluetooth®		•

We reserve the right to make modifications.

seca north america | west
13601 Benson Avenue
Chino, CA 91710
USA

seca north america | east
1439 Stuart Engals Blvd, Suite 200
Mount Pleasant, SC 29464
USA

T +1 800 542 7322
F +1 888 705 7397
E info.us@seca.com

Precision for health
seca.com

12/2024