seca **956**



- The simple solution for disabled people.
- Fold-up armrests and swiveling footrests make sitting down easy.
- Brakes on wheels for secure stance.
- Intuitive operating elements simplify use.
 - All essential functions for maximum usage.



seca 956:

Relief for patient and medical personnel.

Standing on a scale is a problem for physically disabled or older people.

Chair scales are of great assistance to patients and medical staff in nursing homes, hospitals, geriatric and orthopedic departments. The new seca 956 offers a sound, user-oriented solution with all the functions that a good chair scale should offer.



With a capacity of 200 kg and a graduation of 100 grams, the seca 956 fulfills all basic medical standards. Added to that are the typical seca advantages such as fold-up armrests and swiveling footrests, secure wheel brakes, an ergonomically shaped seat and smoothrunning rollers, which, together with the ergonomically shaped handle, make it very simple to push the chair scale.

Simple operation.

Thanks to the fold-up armrests and footrests, elderly and heavy patients can take a seat on the scale easily and comfortably. It is just as easy for personnel to use the scale because the batteries eliminate the need for connection to an electrical outlet and keep on working through up to 4.000 weighings. No time is lost to a learning curve. A user can operate the scale intuitively. The display elements are conveniently located in the handle.



TARE and HOLD.

Of course the seca 956 is also equipped with the important TARE function which allows every newly added weight to be ascertained separately and the HOLD function which lets medical personnel take care of the patient first and read out the results later. All these advantages plus well-known workmanship and durability make the scale a genuine seca product. It fulfills all the expectations people have of the world leader in medical weighing and measuring.







To make sitting down and standing up more convenient, the armrests and footrests fold up out of the way.



The practically positioned display is easy to use.

seca **956**

Technical Data

- Capacity: 200 kg
- Graduation: 100 g
- Dimensions (WxHxD):
 563 x 897 x 956 mm
- Weight: 22.4 kg
- Power supply: batteries
- Functions: TARE, HOLD, automatic switch-off
- Approval class: