

## A range of solutions to optimize patient care in your medical home.

With over 20 million patients being treated at community health centers across the nation, it is essential that every facility is equipped with innovative measuring and weighing devices that support comprehensive preventative care. seca is dedicated to helping optimize and expand quality care to your communities by offering a range of solutions to suit your health center.



*The weighing tray can  
be removed quickly  
and easily at the  
touch of a button.*



### seca 354

Digital baby scale also converts  
to a flat scale for children

Two scales in one: the seca 354 baby scales with a comfortable weighing tray and a fine graduation of 0.5 oz. In no time at all, it can be converted into a floor scale for children with a capacity of 44 lbs. The weighing tray and the base which serves as the floor scale are securely locked together, thus ensuring that the baby can be weighed absolutely safe. A further advantage of the scale is its Breast-Milk-Intake-Function which determines how much milk the baby has consumed. The robust but lightweight scale is easy to transport and self-explanatory buttons make it easy to operate.

- Capacity: 44 lbs / 20 kg
- Graduation: 0.5 oz < 22 lbs > 1 oz / 10 g < 10 kg > 20 g
- Functions: TARE, BMIF, auto-HOLD, lbs/kg switch-over, automatic switch-off
- Optional: carrying case seca 414 and seca 413

# seca | 360° wireless

## One step ahead of the future: Direct EMR integration.

Advance your patient centered medical home by offering a range of solutions designed to embrace health care demands and the increased challenges of your growing community. Only seca's EMR-integrated products like the seca 360° wireless system transmit data securely and help with eliminating manual entry error. Within seconds, measurements are transmitted wirelessly and assigned to the correct patient ID in your EMR system.



### seca 374

Digital baby scale  
with wireless transmission

Ergonomic, convenient and time-saving – those are the words that best describe the seca 374. Its extra large tray lets you weigh a seated toddler of up to 44 lbs, even if the child is fidgety. The adjustable damping function combined with the auto-HOLD function ensures fast and precise work. Net weight of a diapered baby? Choose TARE. The amount of milk breast-fed to a baby? Use the BMIF (Breast-Milk-Intake-Function). Weigh and measure in just one step? Yes – with the optional measuring rod seca 233. Transmit the weight measurement wirelessly to the seca 360° wireless network? Just press SEND. And you can do it wherever you like because the seca 374 is operated by batteries, so there's no need for electrical outlets.

- Capacity: 44 lbs / 20 kg
- Graduation: 0.2 oz < 22 lbs > 0.5 oz / 5 g < 10 kg > 10 g
- Functions: TARE, BMIF, auto-HOLD, auto-CLEAR, lbs/kg switch-over, damping, automatic switch-off, SEND/PRINT
- Optional: measuring rod seca 233, head and foot positioners seca 418, power adapter seca 447, switch-mode power adapter seca 400



The measured height is shown in the display on the headpiece.

### new

### seca 284

Digital measuring station for height  
and weight with wireless transmission

The 360° wireless measuring station seca 284 weighs and measures in a single step. Other benefits include the high-capacity scale integrated in a stable glass platform, the heel positioner and the solid headpiece with the integrated seca Frankfurt Line – for precise head positioning. Not to mention the two displays. One is for a direct read-out of height from the headpiece and the other, a three-line multi-function touch display, shows weight, height and the automatically calculated BMI. On top of all that, the scale shows the slightest weight changes with its fine 50-gram graduation and transmits all measurements via the seca 360° wireless.

#### Scales:

- Capacity: 660 lbs / 300 kg
- Graduation: 0.1 lbs / 50 g
- Functions: TARE, pre-TARE, mother/child-function, HOLD, auto-HOLD, auto-BMI, CLEAR, lbs/kg switch-over, damping, SEND/PRINT

#### Measuring rod:

- Measuring range: 11 inch – 7 ft 2" / 30 – 220 cm
- Graduation: 0.05" / 1 mm
- Functions headpiece: HOLD, user-defined zero setting, automatic switch-off

### optimized



seca 220 measuring rod included.

### seca 703

High capacity column  
scale with wireless transmission

Weighing small children to obese patients – the seca 703 can do it all. Especially when it comes to finding out BMI and height. Then there are two options: integrate a mechanical measuring rod or receive data wirelessly from a separate digital stadiometer. Other advantages include the high load bearing capacity of up to 660 lbs for heavyweight patients, the large and flat platform and time-saving functions Pre-TARE, HOLD, Auto-HOLD and mother/child-function. With everything primed for the future and integration in EMR using the seca 360° wireless.

- Capacity: 660 lbs / 300 kg
- Graduation: 0.1 lbs < 330 lbs > 0.2 lbs / 50 g < 150 kg > 100 g
- Functions: TARE, pre-TARE, mother/child-function, HOLD, auto-HOLD, BMI, auto-CLEAR, damping, automatic switch-off, SEND/PRINT
- Optional: measuring rod seca 224, switch-mode power adapter seca 400

# 360°



### optimized



The scale can be folded for easy transport.

### seca 676

Digital wheelchair scale with handrail and wire-  
less transmission

Stable, functional and mobile: thanks to its large platform and sturdy handrail, the seca 676 is extremely versatile. Patients can be weighed while sitting in a wheelchair or on a chair. The handrail provides valuable support for people who are frail or can walk only with great difficulty. Its high load capacity makes the seca 676 ideal for weighing very heavy patients. After use the scale can be folded together to save space. The sturdy locking device between rail and platform ensures that the seca 676 stands safely even when folded. The handrail also serves as handle when the scale is folded, allowing the scale to be moved around or stowed away effortlessly and quickly on its transport castors. Weighing is simple and user-friendly with the clearly designed, operable display at hip level, additional functions of TARE, auto-HOLD, BMI and the SEND key to transmit measurements wirelessly to a seca wireless printer or the PC.



Integrated ramp folds up for troublefree transport.

- Capacity: 800 lbs / 360 kg
- Graduation: 0.1 lbs / 50 g
- Functions: TARE, pre-TARE, HOLD, auto-HOLD, BMI, auto-CLEAR, lbs/kg switch-over, damping, SEND/PRINT

### Contact us

Toll free +1 800 542 7322 or  
Email: [info.us@seca.com](mailto:info.us@seca.com)

Learn more about  
seca 360° wireless products at  
[www.seca.com/us\\_chc](http://www.seca.com/us_chc)

## seca 869

Flat scale with cable remote display



A two-meter-long cable connects the display with the weighing base. Advantage: the display can be set up separately from the base and read from its position on the wall or desk. Additional benefits include the many extra weighing functions: BMI for nutritional conditions, HOLD to keep measurements displayed and TARE for weighing small children held in an adult's arms. It is suitable for stationary or mobile use, thanks to its handle, compact size, low weight and a leveling base for set-up on all types of surfaces.

- Capacity: 550 lbs / 250 kg
- Graduation: 0.2 lbs < 330 lbs > 0.5 lbs / 100 g < 150 kg > 200 g
- Functions: TARE, HOLD, BMI, lbs/kg switch-over, automatic switch-off
- Optional: stadiometer seca 217 combined with adapter element seca 437, switch-mode power adapter seca 400, carrying case seca 414 (for combination with stadiometer and adapter element), carrying case seca 421

## seca 952

Chair scale for weighing while seated

With a capacity of 440 lbs and a graduation of 0.2 lbs, the seca 952 fulfills all basic medical standards. Added to that are the typical seca advantages such as fold-up armrests and footrests, secure wheel brakes, an ergonomically shaped seat. Its smooth-running rollers and ergonomically shaped handle make it very easy to push the chair scale. Of course the seca 952 is also equipped with the important TARE function which allows every newly added weight to be ascertained separately and the auto-HOLD function which lets medical personnel take care of the patient first and read out the results later.

- Capacity: 440 lbs / 200 kg
- Graduation: 0.2 lbs / 100 g
- Functions: TARE, auto-HOLD, lbs/kg switch-over, automatic switch-off



## seca 216

Mechanical measuring rod for children and adults

The seca 216 has a variable scale with a range of 1.4 to 90 inch that is inserted into a 55.3 inch frame. The plastic frame is simply mounted on the wall at the desired height. As the user can choose the range shown in the frame, the seca 216 is suitable for measuring children or adults. The results are clearly displayed in a window on the side of the rod and held in place by the locking screw on the head piece. Suitable for international use, the scale is printed in centimeters on one side of the tape and centimeters and inches on the other.

- Measuring range: max. 55.3" / 138 cm, variable from 1.4–90" / 3.5–230 cm (freely selectable)
- Graduation: 1/16" / 1 mm

