



# WEIGHING& MEASURING REDEFINED

seca SCALE-UP LINE ———



Play video to see the scale in action!

## Built for Dunability Designed for Safety Ready for High Traffic

MOBILITY

MEETS

DURABILITY

DURABILITY

The optimized mobile and easy-to-move model of the seca Scale-up Line incorporates a digital height rod and optional rechargeable battery pack. Digital height measuring results appear in the display automatically with millimeter precision.

The seca Scale-up Line is purpose-built to endure the rigors of daily medical work with resilience. Its heavy-duty, robust design ensures consistent accuracy even after years of relentless service, providing healthcare teams with trusted reliability.







Obtain weight and digital height conveniently at each patient room, bedside, or wherever



The rock-solid steel handrail enhances stability to reduce fall risk.

### DESIGNED TO THRIVE

#### in Demanding Healthcare Environments

- (1) Ultrasound height measuring
- (2) Sturdy aluminum column
- (3) Customizable touch-display
- (3) Voice-guided automation
- (4) Robust steel handrail
- 5) Patient and user identification\*
- 6) Optional Rechargabel battery pack
- (7) Heavy-duty wheels
- (8) Break-resistant Bearclaw glass
- (8) Easy to clean surfaces
- (+) Capacity of 800lb / 360kg
- + Fine 0.1oz / 50 g graduation



The seca Scale-up Line redefines medical scales' standards, delivering precision and versatility that will reshape your expectations. Its innovative features and unmatched durability set a new benchmark in the industry, making it the go-to choice for healthcare professionals.

#### One Scale to Master all Weighing Challenges

The seca Scale-up Line is your versatile solution, perfect for a wide range of patients and healthcare scenarios. Whether you're working with pediatric or geriatric patients, need mobile functionality, or require seamless EMR integration, this scale was designed to excel in every environment.



Bedside support for convenient weighing

Optional barcode scanner mount seca 463, and barcode scanner from a thirdparty vendor.



Seamless Direct EMR-validated Integration

Enhance efficiency, eliminate transcription errors, and foster a zero-incident environment with the seca Scale-up Line.
Seamlessly integrate patient measurements into your EMR\*
System. Imagine the relief for staff as patient weight and height data is automatically captured.
Contact us for more information on our integration solutions.

→ seca.com/integration










\* EMR integration features are functional exclusively when the device is connected to the AC Adapter and are not applicable when powered by a mobile battery pack.





## TAILOR YOUR SOLUTION

Find the perfect Scale-up Line model for your needs



13601 Benson Avenue Chino, CA 91710 USA

#### seca north america | east

1439 Stuart Engals Blvd, Suite 200 Mount Pleasant, SC 29464 USA

- T +1 800 542 7322
- F +1 888 705 7397
- E info.us@seca.com

seca.com/scaleup

